Urea 37


o-Phthalaldehyde 37°C. Colorimetric

Quantitative determination of urea IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Urea in the sample reacts with o-phthalaldehyde in acid medium forming a coloured complex that can be measured by spectrophotometry:

 The intensity of the color formed is proportional to the urea concentration in the sample¹.

CLINICAL SIGNIFICANCE

Urea is the final result of the metabolism of proteins; it is formed in the liver from their destruction.

 It can appear the urea elevated in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	o-Phthalaldehyde	4,8 mmol/L
R 2	Borate solution Sulphuric acid	87 mmol/L 3 mol/L
UREA CAL	Urea aqueous primary standard 50 mg/dL	

PRECAUTION

R2: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 510 nm $\geq 0,20$.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 510 nm.
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment^(Note 2)

SAMPLES

 - Serum or heparinized plasma¹: Do not use ammonium salts or fluoride as anticoagulants.

 - Urine¹: Dilute sample 1/50 in distilled water. Mix. Multiply the results by 50 (dilution factor). Preserve urine samples at pH < 4.

Urea is stable at 2-8°C for 5 days;

PROCEDURE AND CALCULATIONS

- Assay conditions:
Wavelength: 510 (500-550) nm
Cuvette: 1 cm. light path
Temperature: 37°C
- Adjust the instrument to zero with distilled water.

A) Kinetic method

- Pipette into a cuvette:

	Blank	Standard	Sample
R 1 (mL)	1,0	1,0	1,0
Standard ^(Note 1,3,4) (μL)	--	50	--
Sample (μL)	--	--	50

- Mix, wait 1 minute and add:

R 2 (mL)	1,0	1,0	1,0
----------	-----	-----	-----

- Mix, incubate at 37°C and read the absorbance after 1 minute (A₁) and after 2 minutes (A₂).

- Calculate the increase of the absorbance $\Delta A = A_2 - A_1$.

Calculations

$$\frac{(A_2 - A_1) \text{ Sample} - (A_2 - A_1) \text{ Blank}}{(A_2 - A_1) \text{ Standard} - (A_2 - A_1) \text{ Blank}} \times 50 (\text{Std. conc.}) = \text{mg/dL urea in the sample}$$

B) End point

- Pipette into a cuvette:

	Blank	Standard	Sample
R 1 (mL)	1,0	1,0	1,0
Standard ^(Note 1,3,4) (μL)	--	25	--
Sample (μL)	--	--	25

- Mix, and add:

R 2 (mL)	1,0	1,0	1,0
----------	-----	-----	-----

- Mix and incubate 15 min at 37°C.

- Read the absorbance (A) against the Blank.

Calculations

$$\frac{(A) \text{ Sample} - (A) \text{ Blank}}{(A) \text{ Standard} - (A) \text{ Blank}} \times 50 (\text{Std. conc.}) = \text{mg/dL urea in the sample}$$

$$10 \text{ mg/L urea BUN divided by } 0,466 = 21 \text{ mg/L urea} = 0,36 \text{ mmol/L urea}^1.$$

Conversion factor: mg/dL x 0,1665 = mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum : 15- 45 mg/dL (2,49-7,49 mmol/L)

Urine : 20 - 35 gr/24 h.

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,70 mg/dL to linearity limit of 200 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	43,2	145	41,9	147
SD	1,51	1,10	0,80	2,83
CV (%)	3,49	0,76	1,92	1,91

Sensitivity: 1 mg/dL = 0,00459 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,9884

 Regression equation: $y = 0,975x + 0,6$

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

 It is recommended to use heparin as anticoagulant. Do not use ammonium salts or fluoride¹.

 A list of drugs and other interfering substances with urea determination has been reported by Young et. al^{2,3}.

NOTES

- UREA CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- Glassware and distilled water must be free of ammonia and ammonium salts¹.
- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- Kaplan A. Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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