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**BIBLIOGRAFIA**

1. Clinical and Laboratory Standards Institute (CLSI) H47: One-stage Prothrombin Time test and Activated partial Thromboplastin Time test; approved guideline
2. Clinical and Laboratory Standards Institute (CLSI) H21: Collection, Transport, and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; approved guideline
3. Poller L. The Prothrombin Time. WHO/LAB/98.3. 1998.

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Ref. 1709222		4 x 4 mL
Ref. 1709223	Cont.	10 x 2 mL
Ref. 1709224		10 x 4 mL
Ref. 1709225		6 x 4 mL

**Prothrombin time reagent (PT)**
**IVD**

Store at 2-8°C

**INTENDED USE**

High-sensitivity calcium thromboplastin for determination of Prothrombin Time, for evaluation of the extrinsic coagulation pathway and monitoring Oral Anticoagulant Therapy in human citrated plasma.

**PRINCIPLE**

 Based on A. Quick's test, the assay measures the time elapsed from the mixture of plasma with Thromboplastin (a tissue extract rich in Tissue Factor, phospholipids and calcium) until clot formation. Coagulation is initiated by activation of FVII by Tissue Factor<sup>1</sup>.

**COMPOSITION**

PT	Rabbit brain thromboplastin, calcium chloride, heparin inhibitor and preservatives. Freeze-dried.
----	---

**PRECAUTIONS AND WARNINGS:**

- Avoid contact with skin and eyes. Wear suitable protective clothing.
- The ISI value indicated for each lot of reagent is generic for photo-optical instruments. Local ISI calibration is highly recommended<sup>1</sup>.
- Each laboratory should establish the reference value for each lot of reagent.

**PREPARATION**

Dissolve (→) the contents of each vial of thromboplastin with 4,0 mL (Ref. 1709222-1709224-1709225) or 2 mL (Ref. 1709223) of purified water. Replace the stopper and mix gently. Avoid foaming. Let the vial stand undisturbed for 15 min at room temperature (aprox. 25°C). Mix gently before each use.

**STABILITY**

Unopened reagent is stable until the expiration date shown on the vial, when stored at 2-8°C. Stability after reconstitution: 30 days at 2-8° or 24 hours at 25°C. Do not freeze redissolved product.

**SAMPLES**

 Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Immediately after, centrifuge at 2500 x g, 15 min, and transfer the plasma into siliconized glass or plastic tubes. The samples are stable 24 hours at room temperature (15-25°C). Do not store the samples at 2-8°C to avoid cold-activation of FVII<sup>2</sup>.

**PROCEDURE**

The reagent may be used with any method for end-point (clot) detection, such as manual, photo-optical or mechanical methods/instruments.

**PT test procedure**

1. Prewarm to 37±1°C the amount of Thromboplastin needed and the plasma samples during no more than 10 min.
2. Start the PT test by mixing two parts (example: 200µl) of prewarmed thromboplastin reagent and one part (example: 100µl) of prewarmed citrated plasma. Start a timing device the instant the reagents are mixed. Record the time from final reagent addition to clot formation.

**RESULTS**

Results may be given in seconds, PT rate (comparing with a reference value), percent of activity (%; with calibration curve using a reference plasma (fresh pool or freeze-dried)), or as International Normalized Ratio (INR):

$$INR = \left( \frac{PT (sec)}{PT referenc. (sec)} \right)^{ISI}$$

The ISI (International Sensitivity Index) is lot-specific and it depends on the clot-detection system used. Spinreact indicates for each lot a generic ISI value for photo-optical detection. Such value is traceable to the ISI of the reference thromboplastin RBT/05.

**QUALITY CONTROL**

It is convenient to test control plasmas together with the samples. If the values obtained are outside the tolerance range, the reagents and the clot-detection system should be checked.

COAGULATION CALIBRATOR PLASMA	REF: 1709101
NORMAL CONTROL PLASMA	REF: 1709104
PATHOLOGIC CONTROL PLASMA	REF: 1709106

Each laboratory should establish its own Quality Control program, including correction measures.

**REFERENCE VALUES**

Using the BioBas10 coagulometer and with 20 normal samples the following normal range (95% confidence) was established: 11,1 – 14,3 seconds. Each laboratory should establish its own reference values.

**PERFORMANCE CHARACTERISTICS**
**Precision:** Precision was studied with Spinreact's BioBas1000 using lyophilized normal (level 1) and pathologic (level 2) plasma pools, during 19 days (2 runs/day) and four replicates in each run:

Level	n	Mean (secs.)	Coefficient of Variation		
			Within-Run	Between-Run	Total
1	112	13,8	1,6 %	0 %	2,2 %
2	112	37,3	1,7 %	0 %	2,5 %

**Sensitivity:**

FII		FV		FVII		FX	
%	PT (s)	%	PT (s)	%	PT (s)	%	PT (s)
91	13,3	89	13,2	97	13,2	94	13,5
73	13,5	71	13,8	78	13,4	75	14,1
50	14,1	49	14,7	48	14,2	52	14,8
36	14,6	40	15,1	44	14,4	42	15,8
32	14,9	31	16,0	34	14,7	33	16,5

**INTERFERENCES**

The reagent contains a specific heparin inhibitor. Concentrations ≤1UI/mL do not interfere with the results.

Discard hemolytic or coagulated samples. With coagulometer BioBas1000, triglycerides (Intralipid) don't interfere up to 500 mg/dL, hemoglobin does not interfere up to 100 mg/dL, and bilirubin does not interfere up to 15 mg/dL.

 Some common drugs may affect the PT result <sup>3</sup>.

**BIBLIOGRAPHY**

1. Clinical and Laboratory Standards Institute (CLSI) H47: One-stage Prothrombin Time test and Activated partial Thromboplastin Time test; approved guideline
2. Clinical and Laboratory Standards Institute (CLSI) H21: Collection, Transport, and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; approved guideline
3. Poller L. The Prothrombin Time. WHO/LAB/98.3. 1998.

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