Quantitative determination of apolipoprotein A-II (APO A-II) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Turbidimetric test for the measurement of apolipoprotein A-II in human serum or plasma.

Anti- Apo A-II antibodies when mixed with samples containing Apo A-II, form insoluble complexes. These complexes cause an absorbance change, dependent upon the Apo A-II concentration of the patient sample, that can be quantified by comparison from a calibrator of know Apo A-II concentration.

CLINICAL SIGNIFICANCE¹

Lipoprotein consists of lipid and protein portions, and the protein portion is referred to as apolipoprotein.

Apolipoprotein A-II is known as the main structural protein of HDL (High Density Lipoprotein). In some cases, behavior of apolipoprotein A-II is different from that of apolipoprotein A-I, which is another main structural protein of HDL. Therefore, the ratio of A-I/A-II is drawing increasing attention.

REAGENTS

Diluent (R1)	2-Amino-2-hydroxymethyl-1,3-propanediol buffer 100 mmol/L, pH 8.5. Macrogol
Antibody (R2)	Anti-human apolipoprotein A-II goat-polyclonal antibody
Optional	APO CAL ref: 93005

CALIBRATION

The assay and the value of the calibrator concentration have been standardized against an Internal Reference Material. It is recommended the use of the APO CAL Calibrator for calibration.

PREPARATION

Reagents: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 600 nm filter.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 2 weeks at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolized or lipemic samples.

PROCEDURE

- Bring the reagents and the photometer (cuvette holder) to 37°C.
- Assay conditions:

Wavelength : 600 nm
 Temperature : 37 °C
 Cuvette lighth path : 1cm

- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

Reagent R1 (µL)	300
Sample or Calibrator (µL)	5

- Mix and read the absorbance (A₁) exactly 5 minutes after the sample addition.
- Immediately, pipette into de cuvette:

Reagent R2 (µL)	100
-----------------	-----

- Mix and read the absorbance (A₂) of calibrators and sample exactly 5 minutes after the R2 addition.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

CALCULATIONS

$$\frac{(A_2 - A_1)_{\text{sample}}}{(A_2 - A_1)_{\text{calibrator}}} \times \text{Calibrator concentration} = \text{mg/dL Apo A-II}$$

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact Apolipoprotein Control (Ref.:93006) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES⁵

Between 25.1 – 34.5 mg/dL.

Each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

1. Linearity: Up to 60 mg/dL (Nota 1), under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Detection Limit: Values less than 0.15 mg/dL give non-reproducible results.

3. Sensitivity: Apo A-II Calibrator 31.5 mg/dL = 5.94 ΔmA / mg/dL apo A-II.

4. Precision:

Mean (mg/dL)	Intra-assay (n=10)			Inter-assay (n=5)	
	26.7	40.0	53.7	35.2	43.3
SD	0.6	0.27	0.21	0.38	0.32
CV	2.40	0.67	0.39	1.1	0.7

5. Accuracy: Results obtained using this reagent (y) were compared to those obtained with single immunodiffusion method. 50 samples of Apo A-II were assayed. The correlation coefficient (r) was 0.971 and the regression equation $y = 0.97x + 0.56$ for serum, and $r = 0.978$ $y = 1.02x - 0.22$ for plasma.

The results of the performance characteristics depend on the used analyzer.

INTERFERENCES

Hemoglobin (up to 1000 mg/L), bilirubin (up to 40 mg/dL), and lipemia (up to 5%), do not interfere. Other substances may interfere^{6,7}.

NOTES

- Linearity depends on the calibrator concentration.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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