Quantitative determination of Plasma Proteins
IVD

Store at 2-8°C.

PRODUCT CHARACTERISTICS

The PROT CONTROL is a liquid serum from human origin containing protein concentrations suitable to evaluate the accuracy and precision of plasma protein determinations by immunoturbidimetry method.

REAGENTS

PROT-CONTROL	Control. Human serum. Sodium azide 0.95 g/L. Protein concentrations are stated on the list below.
--------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control are traceable to Standard Reference Material CRM 470/RPHHA (Institute for Reference Materials and Measurements, IRMM) (all constituents except antithrombin III).

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Opened vials, stable for 30 days at 2-8°C .

PROCEDURE

To be used in immunoturbidimetric assays.

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LOT

This Control data sheet for lot is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number. This product was manufactured, packaged and tested in accordance to SPINREACT, S.A./S.A.U. Quality System and meets all product specifications. Our Quality System has been approved by an accredited registering body to ISO 9001.

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Controlo Geral de Proteínas/ General Protein Control

Immunoturbidimetry / Imunoturbidimetria

Quantitative determination of Plasma Proteins IVD

Store at 2-8°C.

PRODUCT CHARACTERISTICS

The PROT CONTROL is a liquid serum from human origin containing protein concentrations suitable to evaluate the accuracy and precision of plasma protein determinations by immunoturbidimetry method.

REAGENTS

PROT-CONTROL	Control. Human serum. Sodium azide 0.95 g/L. Concentrations of the constituents are stated on the list below.
--------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control are traceable to Standard Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM) (all constituents except antithrombin III).

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Opened vials, stable for 30 days at 2-8°C.

PROCEDURE

To be used in immunoturbidimetric assays.

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