
Quantitative determination of α -amylase IVD.

Store at 2-8°C

PRINCIPLE OF THE METHOD

α -Amylase hydrolyzes the 2-chloro-4-nitrophenyl- α -D-maltotriose (CNPG₃) to release 2-chloro-4-nitrophenol (CNP) and form 2-chloro-4-nitrophenyl- α -D-maltoside (CNPG₂), maltotriose (G₃) y glucose (G) according to the following reaction:

The rate of 2-chloro-4-nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of α -amylase present in the sample¹.

CLINICAL SIGNIFICANCE

α -Amylase (AMS) is an enzyme that helps to digest the glycogen and the starch. It is produced mainly by exocrine pancreas and salivary glands. This determination is made mainly in diagnosis or to control diseases of the pancreas as acute or chronic pancreatitis. It can also reflect biliary or gastrointestinal disease and other upheavals^{2,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	MES pH 6.0	100 mmol/L
	CNPG ₃	2.25 mmol/L
	Sodium clorhidre	350 mmol/L
	Calcium acetate	6 mmol/L
	Potassium thiocyanate	900 mmol/L
	Sodium azide	0.95 gr/L
OPTIONAL	Spintról H CAL	

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

After opening, the reagent is stable for 60 days when properly capped immediately after each opening and stored at 2-8°C.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm \geq 0,40.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 405 nm.
- Thermostatic bath at 37°C^(Note 1).
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment^(Note 2).

SAMPLES

- Serum or plasma¹, remove from cells as soon as possible. It is recommended to use heparin as anticoagulant.
 - Urine, adjust pH to approximately 7.0 prior to storage.
- Stability: 1 month at 2-8°C.

REFERENCE VALUES¹

Serum or plasma Up to 90 U/L of α -amylase

Urine Up to 450 U/L of α -amylase

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION SPINLAB 180

Name	Amylase	Ref. male low	5
Abbr. Name	AMI	Ref. male high	90
Mode	Kinetic	Ref. female low	5
Wavelength	405 nm	Ref. female high	90
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	5 U/L	Panic value low	*
High Conc.	2000 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 μ L	normal volume	300 μ L
rerun volume	300 μ L	rerun volume	300 μ L
Sample		Sample	
normal volume	6.0 μ L	normal volume	6.0 μ L
rerun volume	3.0 μ L	rerun volume	3.0 μ L
R2 bottle (mL)	5 mL		
normal volume	0.0 μ L		
rerun volume	0.0 μ L		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	51, 175 sec.
Linearity limit	10.0 %	Linearity limit	10.0 %
Factor		Factor	
Reagent blank	No	Reagent blank	No
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTRÓL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

- α -Amylase enzyme activity is temperature dependent. Assays performed at temperatures $<37^\circ\text{C}$ or $>37^\circ\text{C}$ will show an apparent decrease or increase levels.
- Saliva and sweat contain α -amylase. Avoid mouth pipetting and skin contact with the reagent or material used.
- Contains potassium thiocyanate. Avoid inhalation, skin or eyes contact.
If it happens, wash with plenty of water and consult a doctor.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- Ying Foo A et al. Amylase measurement with 2-chloro-4-nitrophenyl maltotriose as substrate. Clin Chim 272, 1998; 137-147.
- McNeely M. Amylase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1112-116.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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1. Ying Foo A et al. Amylase measurement with 2-chloro-4-nitrophenyl maltotrioxide as substrate. *Clin Chim* 272, 1998; 137-147.
2. McNeely M. Amylase. Kaplan A et al. *Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton* 1984; 1112-116.
3. Young DS. *Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.*
4. Young DS. *Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.*
5. Burtis A et al. *Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.*
6. Tietz N W et al. *Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.*

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 Ref: SP41201 R: 10 x 25 mL