
Quantitative determination of lactate dehydrogenase (LDH)
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Lactate dehydrogenase (LDH) catalyses the reduction of pyruvate by NADH, according the following reaction:

The rate of decrease in concentration of NADPH, measured photometrically, is proportional to the catalytic concentration of LDH present in the sample¹.

CLINICAL SIGNIFICANCE

Lactate dehydrogenase (LDH) is an enzyme with wide tissue distribution in the body.

The higher concentrations of LDH are found in liver, heart, kidney, skeletal muscle and erythrocytes.

Increased levels of the enzyme are found in serum in liver disease, myocardial infarction, renal disease, muscular dystrophy and anemia^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

Reagent 1	Phosphate pH 7.8	80 mmol/L
Buffer	Pyruvate	0,6 mmol/L
Reagent 2	NADH	0,18 mmol/L
Substrate		
Optional	SPINTROL H CAL	

PREPARATION

DUAL MODE: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm <1.00.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C o 37°C (± 0.1°C)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum¹. Separated from cells as rapidly as possible. Do not use oxalates as anticoagulants since they inhibit the enzyme. Do not use haemolysed samples. Stability: 2 days at 2-8°C.

INTERFERENCES

Haemolysis interferes with the assay.

Some anticoagulants such as oxalates interfere with the reaction¹.

A list of drugs and other interfering substances with LDH determination has been reported^{2,3}.

APPLICATION SPINLAB 180

Name	LDH	Ref. male low	230 U/L
Abbr. Name	LDH	Ref. male high	460 U/L
Mode	Kinetic	Ref. female low	230 U/L
Wavelength	340 nm	Ref. female high	460 U/L
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	5 U/L	Panic value low	*
High Conc.	1450 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000

DUAL MODE

Sample blank	No
R1 bottle (mL)	25 mL
normal volume	240 µL
rerun volume	240 µL
Sample	
normal volume	5.0 µL
rerun volume	5.0 µL
R2 bottle (mL)	5 mL
normal volume	60.0 µL
rerun volume	60.0 µL
Predilución	No
Slope blank	No
Delay, min. time	50, 186 sec.
Linearity limit	10.0 %
Factor	
Reagent blank	No
Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

25°C	30°C	37°C
120-240 U/L	160-320 U/L	230-460 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

1. Pesce A. *Lactate dehydrogenase*. Kaplan A et al. *Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton* 1984; 1124-117, 438.
2. Young DS. *Effects of drugs on Clinical Lab. Tests*, 4th ed AACC Press, 1995.
3. Young DS. *Effects of disease on Clinical Lab. Tests*, 4th ed AACC 2001.
4. Burtis A et al. *Tietz Textbook of Clinical Chemistry*, 3rd ed AACC 1999.
5. Tietz N W et al. *Clinical Guide to Laboratory Tests*, 3rd ed AACC 1995.

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