## Quantitative determination of albumin IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Albumin in the presence of bromcresol green at a slightly acid pH, produces a colour change of the indicator from yellow-green to green-blue. The intensity of the color formed is proportional to the albumin concentration in the sample<sup>1,2,3,4</sup>.

### CLINICAL SIGNIFICANCE

One of the most important serum proteins produced in the liver is albumin.

This molecule has an extraordinarily wide range of functions, including nutrition, maintenance of oncotic pressure and transport of Ca<sup>++</sup>, bilirubin, free fatty acid, drugs and steroids. Variation in albumin levels indicate liver diseases, malnutrition, skin lesions such as dermatitis and burns or dehydration<sup>1,7,8</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R</b>	Bromcresol green pH 4.2	0,12 mmol/L
<b>Optional</b>	SPINTROL H CAL	

### PREPARATION

Ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 630 nm  $\geq$  0.40.

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 630 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

### SAMPLES

Serum or plasma, free of hemolysis<sup>1</sup>: Stability 1 month at 2-8°C or 1 week at 15-25°C.

### REFERENCE VALUES

3.5 to 5.0 g/dL<sup>1</sup>.

These values are for orientation purpose; each laboratory should establish its own reference range.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## APPLICATION SPINLAB 180

Name	Albumin	Ref. male low	3,5
Abbr. Name	ALB	Ref. male high	5,0
Mode	Endpoint	Ref. female low	3,5
Wavelength	620 nm	Ref. female high	5,0
Units	g/dL	Ref. Ped. Low	3,8
Decimals	1	Ref. Ped. High	5,5
Low Conc.	0,0 g/dL	Control 1	*
High Conc.	6,0 g/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	398 $\mu$ L	normal volume	398 $\mu$ L
rerun volume	398 $\mu$ L	rerun volume	398 $\mu$ L
Sample		Sample	
normal volume	2,0 $\mu$ L	normal volume	2,0 $\mu$ L
rerun volume	2,0 $\mu$ L	rerun volume	2,0 $\mu$ L
R2 bottle (mL)	5 mL		
normal volume	0 $\mu$ L		
rerun volume	0 $\mu$ L		
Predilución	No		
Incubation time	11.5 min	Incubation time	11.5 min
Factor		Factor	
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 0,0349 g/dL to linearity limit of 6 g/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

Mean (g/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	5,00	3,71	4,56	3,07
SD	0,02	0,02	0,28	0,18
CV (%)	0,47	0,55	6,20	5,90

**Sensitivity:** 1 g/dL = 0,2003 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99169.

Regression equation: y = 1,045x - 0,028.

The results of the performance characteristics depend on the analyzer used.

### NOTES

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

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1. Gendler S. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1268-1273 and 425.
2. Rodkey F L. Clin Chem 1965; 11: 478-487.
3. Webster D. Clin Chem. 1974; Acta 53: 109-115.
4. Doumas BT Clin Chem. 1971; Acta 31: 87-96.
5. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
6. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
7. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
8. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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