Protein in urine and CSF

Pyrogallol red. Colorimetric

Quantitative determination of total urinary and CSF protein.

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Protein react in acid solution with pyrogallol red and molybdate to form a colored complex.

The intensity of the color formed is proportional to the protein concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

In healthy persons, the urine contains no protein or only a trace amount of protein; normally the glomeruli prevent passage of protein from the blood to the glomerular filtrate. Glomerular injury causes increased permeability to plasma proteins, resulting in proteinuria, which refers to the presence of protein in the urine.

A persistent finding of proteinuria is the single most important indication of renal disease.

Elevated concentration of protein in cerebro-spinal fluid (CSF) can be caused by infections and intracranial pressure^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Pyrogallol red	50 µmol/L
	Sodium molybdate	0,04 mmol/L
PROTEIN U & CSF CAL	Albumin/Globulin aqueous primary standard 1000 mg/L	

PREPARATION

The reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 598 nm $\geq 0,30$.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 598 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

- Urine 24 h: Stability 8 days at 2-8°C.
- Cerebrospinal fluid (CSF): Stable 4 days at 2-8°C.

REFERENCE VALUES⁵

Urine:	< 100 mg/24 h (< 150 mg/24 h in pregnancy)
CSF:	Children 300 -1000 mg/L
	Adults 150 - 450 mg/L

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION SPINLAB 180

Name	PROT.ORINA	Ref. male low	**
Abbr. Name	PTOR	Ref. male high	**
Mode	Endpoint	Ref. female low	**
Wavelength	578 nm	Ref. female high	**
Units	mg/L	Ref. Ped. Low	**
Decimals	0	Ref. Ped. High	**
Low Conc.	100 mg/L	Control 1	*
High Conc.	4000 mg/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	10.0 µL	normal volume	10.0 µL
rerun volume	8.0 µL	rerun volume	8.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
incubation time	11.5 min.	Incubation time	11.5 min.
Factor		Factor	
Reagent blank	Si (0.000)	Reagent blank	Si (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

QUALITY CONTROL

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.
3. **Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Orsonneau JL et al. An improved Pyrogallol Red-Molybdate Method for Determining Total Urinary Protein. Clin Chem 1989 (35):2233-2236.
2. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis, Toronto. Princeton 1984; 1316-1324 and 418.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Name	PROT.ORINA	Ref. male low	**
Abbr. Name	PTOR	Ref. male high	**
Mode	Endpoint	Ref. female low	**
Wavelength	578 nm	Ref. female high	**
Units	mg/L	Ref. Ped. Low	**
Decimals	0	Ref. Ped. High	**
Low Conc.	100 mg/L	Control 1	*
High Conc.	4000 mg/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	10.0 µL	normal volume	10.0 µL
rerun volume	8.0 µL	rerun volume	8.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volumen	0 µL		
Predilución	No	Incubation time	11.5 min.
incubation time	11.5 min.		
Factor		Factor	
Reagent blank	Si (0.000)	Reagent blank	Si (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

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1. Orsonneau JL et al. An improved Pyrogallol Red-Molybdate Method for Determining Total Urinary Protein. Clin Chem 1989; 35:2233-2236.
2. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Name	PROT.ORINA	Ref. male low	**
Abbr. Name	PTOR	Ref. male high	**
Mode	Endpoint	Ref. female low	**
Wavelength	578 nm	Ref. female high	**
Units	mg/L	Ref. Ped. Low	**
Decimals	0	Ref. Ped. High	**
Low Conc.	100 mg/L	Control 1	*
High Conc.	4000 mg/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	10.0 µL	normal volume	10.0 µL
rerun volume	8.0 µL	rerun volume	8.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volumen	0 µL		
Predilución	No		
incubation time	11.5 min.	Incubation time	11.5 min.
Factor		Factor	
Reagent blank	Si (0.000)	Reagent blank	Si (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

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