Quantitative determination of phosphorus IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Direct method for determining inorganic phosphate. Inorganic phosphate reacts in acid medium with ammonium molybdate to form a phosphomolybdate complex with yellow color. The intensity of the color formed is proportional to the inorganic phosphorus concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Phosphorus is an essential mineral for tissue bone formation and is required by every cell in the body for normal function. Approximately 85% of the body phosphorus is found in bone and in teeth.

Low levels of phosphorus, can be caused by hypervitaminosis D, primary hyperparathyroidism, renal tubular disorders, antacids or malabsorption.

High levels of phosphorus can be caused by diet, bone metastases, liver disease, alcohol ingestion, diarrhea and vomiting^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Ammonium molybdate	0.40 mM
Molybdc	Sulphuric acid (SO ₄ H ₂)	210 mM
	Detergents	
Optional	Spintrol H CAL	

PRECAUTIONS

H314-Causes severe skin burns and eye damage. Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm \geq 0.54.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment ^(Note 1).

NOTES

- Most of the detergents and water softening products used in the laboratories contain chelating agents and phosphates. It is recommended to rinse glassware in diluted nitric acid and water before using.
- Use clean disposable pipette tips for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

SAMPLES

- Serum or plasma^{1,5}:
Free of hemolysis. Serum or plasma should be removed from the clot as quickly as possible to avoid elevation of serum phosphorus from hydrolysis or leakage of phosphate present in erythrocytes. Stability: 7 days at 2-8°C.
- Urine^{1,2} (24 h):
Collect the specimen into a bottle containing 10 mL of 10% v/v hydrochloric acid (HCl) to avoid phosphate precipitations. Adjust to pH 2. Dilute the sample 1/10 with distilled water. Mix. Multiply the result by 10 (dilution factor). Stability: 10 days at 2-8°C.

APLICACION SPINLAB 180

Name	Phosphorus	Ref. male low	2.5 mg/dL
Abbr. Name	FOSF	Ref. male high	5.0 mg/dL
Mode	Endpoint	Ref. female low	2.5 mg/dL
Wavelength	340 nm	Ref. female high	5.0 mg/dL
Units	mg/dL	Ref. Ped. Low	4.0 mg/dL
Decimals	1	Ref. Ped. High	7.0 mg/dL
Low Conc.	0.3 mg/dL	Control 1	*
High Conc.	15 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat.offset	0.000

DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3.0 µL	normal volume	3.0µL
rerun volume	2.0 µL	rerun volume	2.0µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
incubation time	4.5 min.	Incubation time	4.5 min.
Factor		Factor	
Reagent blank	Si (0.000)	Reagent blank	Si (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum or plasma:

Children	4,0 – 7,0 mg/dL	≅	1,29 – 2.26 mmol/L
Adults	2,5 – 5,0 mg/dL	≅	0.80 – 1.61 mmol/L

Urine:

Adults	0.4 – 1.3 g /24 h
--------	-------------------

These values are for orientation purpose; each laboratory should establish its own reference range.

BIBLIOGRAPHY

- Farrell E C. Phosphorus. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1072-1074 and 418.
- Daly J A. et al. Clin Chem 1972; 18 (3): 263-265.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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Name	Fosforo	Ref. male low	2.5 mg/dL
Abbr. Name	FOSF	Ref. male high	5.0 mg/dL
Mode	Endpoint	Ref. female low	2.5 mg/dL
Wavelength	340 nm	Ref. female high	5.0 mg/dL
Units	mg/dL	Ref. Ped. Low	4.0 mg/dL
Decimals	1	Ref. Ped. High	7.0 mg/dL
Low Conc.	0.3 mg/dL	Control 1	*
High Conc.	15 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3.0 µL	normal volume	3.0 µL
rerun volume	2.0 µL	rerun volume	2.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
incubation time	4.5 min.	Incubation time	4.5 min.
Factor		Factor	
Reagent blank	Si (0.000)	Reagent blank	Si (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

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1. Farrell E C. Phosphorus. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1072-1074 and 418.
2. Daly J A. et al. Clin Chem 1972; 18 (3): 263-265.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
5. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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