


**Quantitative determination of cholesterol IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The cholesterol present in the sample originates a coloured complex, according to the following reactions:

The intensity of the color formed is proportional to the cholesterol concentration in the sample<sup>1,2</sup>.

**CLINICAL SIGNIFICANCE**

Cholesterol is a fat-like substance called a lipid that is found in all body cells. The liver makes all of the cholesterol the body needs to form cell membranes and to make certain hormones. The determination of serum cholesterol is one of the important tools in the diagnosis and classification of lipemia. High blood cholesterol is one of the major risk factors for heart disease<sup>5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	PIPES pH 6.9	90 mmol/L
	Phenol	26 mmol/L
	Cholesterol esterase (CHE)	1000 U/L
	Cholesterol oxidase (CHOD)	300 U/L
	Peroxidase (POD)	650 U/L
	4 - Aminophenazone (4-AP)	0.4 mmol/L
<b>Optional</b>	SPINTROL H CAL	

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm ≥ 0.26.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or colorimeter measuring at 505 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

**SAMPLES**

Serum or plasma<sup>1,2</sup>. Stability of the sample 7 days at 2-8°C or freezing at -20°C will keep samples stable for 3 months.

**REFERENCE VALUES**

Risk evaluation<sup>5,6</sup>:

Less than 200 mg/DL	Normal
200-239 mg/dL	Borderline
240 mg/dL and above	High

These values are for orientation purpose; each laboratory should establish its own reference range.

**APPLICATION SPINLAB 180**

Name	Cholesterol	Ref. male low	120.0 mg/dL
Abbr. Name	CHOL	Ref. male high	220.0 mg/dL
Mode	Endpoint	Ref. female low	120.0 mg/dL
Wavelength	505 nm	Ref. female high	220.0 mg/dL
Units	mg/dL	Ref. Ped. Low	0.0 mg/dL
Decimals	1	Ref. Ped. High	220.0 mg/dL
Low Conc.	0.0 mg/dL	Control 1	*
High Conc.	600.0 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3.0 µL	normal volume	3.0 µL
rerun volume	2.0 µL	rerun volume	2.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
Incubation time	4.5 min.	Incubation time	11.5 min.
Factor		Factor	
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**NOTES**

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

**BIBLIOGRAPHY**

1. Naito H.K. *Cholesterol*. Kaplan A et al. *Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.*
2. Meittini F. et al. *The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.*
3. Young DS. *Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.*
4. Young DS. *Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.*
5. Burtis A et al. *Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.*
6. Tietz N W et al. *Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.*

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1. Naito H.K. *Cholesterol*. Kaplan A et al. *Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.*
2. Meiatini F. et al. *The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.*
3. Young DS. *Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.*
4. Young DS. *Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.*
5. Burtis A et al. *Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.*
6. Tietz N W et al. *Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.*

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Ref: SP41021  Cont.  R: 10 x 25 mL