
**Quantitative determination of chloride ion**
**IVD.**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of a red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed is proportional to the chloride ion concentration in the sample<sup>1,2,3,4</sup>.

**CLINICAL SIGNIFICANCE**

It is important clinically the determination of chloride due regulation of osmotic pressure of extra cellular fluid and to its significant role in acid-base balance. Increases in chloride ion concentration may be found in severe dehydration, excessive intake of chloride, severe renal tubular damage and in patients with cystic fibrosis.

Decrease in chloride ion concentration may be found in metabolic acidosis, loss from prolonged vomiting and chronic pyelonephritis<sup>2,7,8</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	Mercuric thiocyanate	2 mmol/L
	Mercuric nitrate	40 mmol/L
	Nitric acid	0.15 mmol/L
		45 mmol/L
<b>Optional</b>	SPINTROL H CAL	

**PRECAUTIONS**

Corrosive (C):R35:Causes severe burns.

Mercury(II) thiocyanate :Harmful (Xn): R20/21/22:Harmful by inhalation, in contact with skin and if swallowed. R33: Danger of cumulative effects.

S13: Keep away from food, drink and animal feeding stuffs. S28: After contact with skin, wash immediately with plenty of water. S45: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). S60: This material and its container must be disposed of as hazardous. S61: Avoid release to the environment. Refer to special instructions/safety data sheets.

**PREPARATION**

Ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

**Do not use reagents over the expiration date.**

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 480 nm  $\geq$  0.15.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or colorimeter measuring at 480 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment (Note 1, 2, 3).

**SAMPLES**

- Serum, plasma, CSF, sweat and other body fluids<sup>1,2</sup>: Free of hemolysis and separated from cells as rapidly as possible. Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.

- Urine<sup>1</sup>: Collect 24-hour urine specimen in chloride free containers. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the sample: Ion chloride is stable 1 week at room temperature (15-25°C), in refrigerator (2-8°C) or frozen (-20°C) temperatures.

**APPLICATION SPINLAB-180**

Name	Chloride	Ref. male low	98
Abbr. Name	CL	Ref. male high	107
Mode	Endpoint	Ref. female low	98
Wavelength	505 nm	Ref. female high	107
Units	mmol/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	0 mmol/L	Control 1	*
High Conc.	130 mmol/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3.0 µL	normal volume	3.0 µL
rerun volume	2.0 µL	rerun volume	2.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
Incubation time	4.5 min.	Incubation time	11.5 min.
Factor		Factor	
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:	95 - 115 mmol/L	CSF:	95 - 110 mmol/L
Urine:	110 - 250 mmol/24h	Sweat:	Up to 60 mmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

**NOTES**

1. It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with H<sub>2</sub>SO<sub>4</sub> - K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> Solution and then thoroughly rinsed it with distilled water.
2. Most of the detergents and water softening products used in the laboratories contains chelating agents. A defective rinsing will invalidate the procedure.
3. Avoid the contact with metal materials.
4. Use clean disposable pipette tips for its dispensation.
5. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

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