Quantitative determination of calcium IVD.

Store at 2-8°C

PRINCIPLE OF THE METHOD

Calcium with Arsenazo III (1,8-Dihydroxy-3,6-disulpho-2,7-naphthalene-bis (azo)-dibenzene arsonic acid), at neutral pH, yields a blue colored complex.

The intensity of the colour formed is proportional to the calcium concentration in the sample^{1,2,3}.

CLINICAL SIGNIFICANCE

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones.

A decrease in albumin level causes a decrease in serum calcium. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption.

Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcosidosis, thyrotoxicosis, hyperparathyroidism^{1,6,7}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Imidazol Buffer pH 6.5	100 mmol/L
Arsenazo III	Arsenazo III	120 mmol/L
Optional	SPINTROL H CAL	

PRECAUTIONS

R: H360- May damage fertility or the unborn child.

CAL: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations are prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 650 nm ≥ 0.50 .

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 650 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment (Note 1, 2).

NOTES

1. It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with diluted 1/1 HNO₃ in water and then thoroughly rinsed it with distilled water.
2. Most of the detergents and water softening products used in the laboratories contains chelating agents. A defective rinsing will invalidate the procedure.
3. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
4. Use clean disposable pipette tips for its dispensation.
5. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

SAMPLES

- Serum or plasma¹: Separated from cells as rapidly as possible. Blood anticoagulants with oxalate or EDTA are not acceptable since these chemicals will strongly chelate calcium.

- Urine¹: Collect 24 hour urine specimen in calcium free containers. The collecting bottles should contain 10 ml of diluted Nitric acid (50% v/v). Record the volume.

Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the samples: Calcium is stable 10 days at 2-8°C.

APPLICATION SPINLAB- 180

Name	Calcium	Ref. male low	9.0
Abbr. Name	CA	Ref. male high	11.0
Mode	Endpoint	Ref. female low	9.0
Wavelength	620 nm	Ref. female high	11.0
Units	mg/dL	Ref. Ped. Low	10.0
Decimals	1	Ref. Ped. High	12.0
Low Conc.	2.0 mg/dL	Control 1	*
High Conc.	30.0 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3.0 µL	normal volume	3.0 µL
rerun volume	2.0 µL	rerun volume	2.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
Incubation time	4.5 min	Incubation time	11.5 min.
Factor		Factor	
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum or plasma:

Adults	8.5-10.5 mg/dL	≅ 2.1-2.6 mmol/L
Children	10 -12 mg/dL	≅ 2.5 - 3 mmol/L
Newborns	8 -13 mg/dl	≅ 2 - 3.25 mmol/L

Urine:

Adults	50 - 300 mg/24h	≅ 1.25 - 7.5 mmol/24h
Children	80 -160 mg/24h	≅ 2 - 4 mmol/24h

These values are for orientation purpose; each laboratory should establish its own reference range.

BIBLIOGRAPHY

1. *Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.*
2. *Kessler G. et al. Clin Chem Vol 10, No 8 1964; 686-706.*
3. *Connerty H. V. et al. Am J Clin Path Vol 45, No 3 1996; 200-296.*
4. *Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.*
5. *Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.*
6. *Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.*
7. *Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.*

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1. *Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.*
2. *Kessler G. et al. Clin Chem 1964; 10 (8); 686-706.*
3. *Connerty H. V. et al. Am J Clin Path 1996; 45 (3); 200-296.*
4. *Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.*
5. *Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.*
6. *Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.*
7. *Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.*

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Ref: SP1001065

Cont.

R: 10 x 25 mL