Quantitative determination of total bilirubin IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Bilirubin (both conjugated and unconjugated) couples with the diazo reagent in the presence of a surfactant to form azobilirubin. The intensity of color formed is proportional to the bilirubin concentration in the sample tested. The increase of absorbance at 546 nm is directly proportional to the total bilirubin concentration.

CLINICAL SIGNIFICANCE

Bilirubin is caused by the degradation of hemoglobin and exists in two forms. Unconjugated bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted. Hyperbilirubinemia is the result of an increase of bilirubin in plasma.

Possible causes:

Total bilirubin: Increase hemolysis, genetic alteration, neonatal anemia, erythropoiesis alterations and presence of drugs.

Direct Bilirubin: cholestasis liver, liver abnormalities and genetic.

Clinical diagnosis should not be made based on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Surfactants Hydrochloric acid (HCl)	<1% 160 mM
R 2	2,4-DPD Hydrochloric acid (HCl) Surfactant	≥2 mM 120 mM <1%

PRECAUTIONS

R1/ R2: H290- Corrosive to metals. H314 - Irritation or skin corrosion.

R1: contains HCl and Triton X-114. R2: contains HCl and 2,4-DPD.

Follow the safety advice given in MSDS and product label.

PREPARATION

The reagents are provided in a ready to use format.

STORAGE AND STABILITY

The reagents are stable until the expiry date stated on the label when stored at 2-8°C, protected from light and contaminations are prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 546 nm.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis. Protect samples from light.

Stability of the sample: 4 days at 2-8°C or 2 months at -20°C.

REFERENCE VALUES

Total bilirubin 0,2-1,2 mg/dL (3,4 – 20,5 μmol/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

Interferences from hemolysis, lipemia and a. ascorbic were evaluated for this total bilirubin method on a Spintech 240 analyzer. Two concentrations of total bilirubin were evaluated. No interferences were observed for lipemia (Intralipid) up to 1800 mg/dL, hemoglobin up to 2000 mg/dL and ascorbic acid up to 40 mg/L.

A list of drugs and other interfering substances with bilirubin has been reported by Young et. al ^{4,5}.

APPLICATION SPINLAB 180

Name	BILIRUBIN T	Ref. male low	*
Abbr.name	BILT	Ref. male high	*
Mode	Twopoint	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/dL	Ref. Ped. low	*
Decimals	2	Ref. Ped. high	*
Low conc.	0.1 mg/dL	Valor pánico bajo	*
High conc.	30.0 mg/dL	Valor pánico alto	*
Calibrator	*	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
Normal volume	240 μL	Normal volume	** μL
Rerun volume	240 μL	Rerun volume	** μL
Sample		Sample	
Normal volume	12.0 μL	Normal volume	** μL
Rerun volume	12.0 μL	Rerun volume	** μL
R2 bottle (mL)	5 mL		
Normal volume	60 μL		
Rerun volume	60 μL		
Predilution	No		
Incubation time	-3.236 seg.	Incubation time	**
Factor	**	Factor	**
Reagent blank	Yes (0,000)	Reagent blank	Yes (0,000)
Low Absorbance	-0.100 Abs Abs	Low Absorbance	-0.100 Abs Abs
High Absorbance	3.000 Abs Abs	High Absorbance	3.000 Abs Abs
R.Abs.L.Limit	-0.100 Abs Abs	R.Abs.L.Limit	-0.100 Abs Abs
R.Abs.H.Limit	3.000 Abs Abs	R.Abs.H.Limit	3.000 Abs Abs

PERFORMANCE CHARACTERISTICS

Measuring range: From *quantification limit* of 0,1 mg/dL to *linearity limit* of 30 mg/dL.

If the results obtained were greater than the linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Inter assay (n= 40)		Intra assay (n= 80)	
		1,169	5,0485	1,1682
SD	0,0285	0,0594	0,012	0,046
CV (%)	2,4	1,2	1,0	0,9

Sensitivity: 1 mg/dL = 0,033 Abs. units.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x) on a Spintech 240 analyzer. The results obtained using 61 samples ranging from 0,42 a 19,36 mg/dL (7,18 to 331,056 μmol/L) were:

Correlation coefficient: (r) 0,996

Regression equation: $y = 0,9836x + 0,1644$

The results of the performance characteristics depend on the analyzer used.

NOTES

1. SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

1. David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307–328.
2. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
3. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP1001046	Cont.	R 1: 10 x 20 mL
		R 2: 10 x 5 mL

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Name	BILIRUBIN T	Ref. male low	*
Abbr.name	BILT	Ref. male high	*
Mode	Twopoint	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/dL	Ref. Ped. low	*
Decimals	2	Ref. Ped. high	*
Low conc.	0.1 mg/dL	Valor pánico bajo	*
High conc.	30.0 mg/dL	Valor pánico alto	*
Calibrator	*	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
Normal volume	240 µL	Normal volume	** µL
Rerun volume	240 µL	Rerun volume	** µL
Sample		Sample	
Normal volume	12.0 µL	Normal volume	** µL
Rerun volume	12.0 µL	Rerun volume	** µL
R2 bottle (mL)	5 mL		
Normal volume	60 µL		
Rerun volume	60 µL		
Predilution	No		
Incubation time	-3.236 seg.	Incubation time	**
Factor	**	Factor	**
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs Abs	Low Absorbance	-0.100 Abs Abs
High Absorbance	3.000 Abs Abs	High Absorbance	3.000 Abs Abs
R.Abs.L.Limit	-0.100 Abs Abs	R.Abs.L.Limit	-0.100 Abs Abs
R.Abs.H.Limit	3.000 Abs Abs	R.Abs.H.Limit	3.000 Abs Abs

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BIBLIOGRAFIA

1. David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
2. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
3. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP1001046

Cont.

R 1: 10 x 20 ml
R 2: 10 x 5 ml