IgA

# Immunoglobulin A

Turbidimetry

## Quantitative determination of human Immunoglobulin A (IgA) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Anti-human IgA antibodies when mixed with samples containing IgA, form insoluble complexes. These complexes cause an absorbance change, dependent upon the IgA concentration of the patient sample, that can be quantified by comparison from a calibrator of known IgA concentration.

### CLINICAL SIGNIFICANCE

IgA represents approximately 10 to 15% of total serum immunoglobulins. Its structure is monomeric, similar to the IgG molecule, but 10 to 15% of IgA in serum is polymeric, particularly IgA<sub>2</sub>, which is more resistant to destruction by some pathogenic bacteria. Another more important form of IgA is called secretory IgA. It is found in tears, sweat, saliva, milk and gastrointestinal and bronchial secretions.

IgA is generally increased in skin, pulmonary, kidney infections, and hepatic cirrhosis. Increased monoclonal IgA concentrations may be found in multiple myeloma and other disturbances of plasmatic cells.

### REAGENTS

<b>R 1</b> Diluent	Tris buffer 20 mmol/L, PEG 8000, pH, 8,3. Sodium azide 0,95 g/L.
<b>R 2</b> Antibody	Goat serum, anti-human IgA, pH 7,5. Sodium azide 0,95 g/L.
<b>Optional</b>	Ref: 1102003 PROT CAL.

### CALIBRATION

The assay is calibrated to the Reference Material CRM 470/RPPHS (Institute of Reference of Materials and Measurements, IRMM). It must be used the PROT CAL Calibrator to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

### PREPARATION

Reagents: Ready to use.

**Calibration Curve:** Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the IgA calibrator by the corresponding factor stated in table below to obtain the IgA concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
ClNa 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1,0

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

### ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

### SPINLAB 180 APPLICATION

Name	IgA	Ref. male low	70 mg/dL
Abbr. Name	IgA	Ref. male high	400 mg/dL
Mode	Twopoints	Ref. female low	70 mg/dL
Wavelength	620 nm	Ref. female high	400 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	1 mg/dL	Control 1	*
High Conc.	600 mg/dL	Control 2	*
Calibrator name	GP CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
<b>Sample blank</b>	<b>No</b>		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	2.0 µL		
rerun volume	2.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	-3.130 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

### REFERENCE VALUES

Between 70 - 400 mg/dL. Each laboratory should establish its own reference range.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

1. *Clinical Guide to Laboratory Tests*, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. *Clin Chem* 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. *Methods in Clinical Chemistry*. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. *Eur J Clin Chem Clin Biochem* 1996; 34: 517-520.
5. Young DS. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Pres, 1997
6. Friedman and Young. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Pres, 1997.

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