Quantitative determination of Immunoglobulin G (IgG) IVD

Store at 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human IgG antibodies when mixed with samples containing IgG, form insoluble complexes. These complexes cause an absorbance change, dependent upon the IgG concentration of the patient sample, that can be quantified by comparison from a calibrator of known IgG concentration.

CLINICAL SIGNIFICANCE

IgG is the most important immunoglobulin produced by plasma cells, and represents about 75% of the total immunoglobulins. Its main function is to neutralize toxins in tissue spaces.

IgG deficit may be due to a congenital primary disturbance (immunodeficiency congenital and acquired) and is a special risk in children.

Polyclonal hyperimmunoglobulinemia is the normal response to infections, especially in hepatitis and cirrhosis as well as autoimmune diseases.

Increases of monoclonal IgG are found in multiple myeloma, lymphocytic leukemia, and Waldenström macroglobulinemia.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Sodium azide 0.95 g/L.
R 2 Antibody	Goat serum, anti-human IgG, pH 7.5. Sodium azide 0.95 g/L.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM). It must be used the PROT CAL Calibrator to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the IgG calibrator by the corresponding factor stated in table below to obtain the IgG concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SPINLAB 180 APPLICATION

Name	IgG	Ref. male low	700 mg/dL
Abbr. Name	IgG	Ref. male high	1600 mg/dL
Mode	Twopoints	Ref. female low	700 mg/dL
Wavelength	620 nm	Ref. female high	1600 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	10 mg/dL	Control 1	*
High Conc.	3000 mg/dL	Control 2	*
Calibrator name	CAL PS	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	2.0 µL		
rerun volume	2.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	-3.130 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 700 - 1600 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. *Clinical Guide to Laboratory Tests*, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. *Clin Chem* 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. *Methods in Clinical Chemistry*. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. *Eur J Clin Chem Clin Biochem* 1966; 14: 401-406.
5. Young DS. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Pres, 1997
6. Friedman and Young. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

1. *Clinical Guide to Laboratory Tests*, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. *Clin Chem* 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. *Methods in Clinical Chemistry*. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. *Eur J Clin Chem Clin Biochem* 1996; 34: 517-520.
5. Young DS. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Pres, 1997
6. Friedman and Young. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Pres, 1997.

PRESENTACION

Ref: SP1103004 Cont. R1: 2 x 20 mL
R2: 2 x 5 mL