Quantitative determination of α₁-antitrypsin (α₁-ATRYP)

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti- α₁-antitrypsin antibodies when mixed with samples containing α₁-antitrypsin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α₁-antitrypsin concentration of the patient sample, that can be quantified by comparison from a calibrator of known α₁-antitrypsin concentration.

CLINICAL SIGNIFICANCE

α₁-antitrypsin is a glycoprotein synthesized in the hepatic parenchyma cells that circulates in the bloodstream. It is the second highest proteinase inhibitor in plasma after α₂-macroglobulin. α₁-antitrypsin is a strong reactor with elastase, skin collagenase, chymotrypsin, plasmin, and thrombin, and also shows inhibitory activity against fungal and leukocytic proteases.

α₁-antitrypsin deficiency is an inherited disorder, and occurs when both parents pass on an abnormal gene (PiZ) to their child. This deficiency is associated with a very high risk for development of pulmonary emphysema and diseases of the liver - neonatal cholestasis, hepatitis, cirrhosis and hepatocellular carcinoma.

α₁-antitrypsin increases in inflammatory or necrosis process. Serum levels begin to rise after approximately 24 hours and peak at 3 or 4 days if the insult is acute and short-lived.

REAGENTS

| | |
|-----------------|---|
| R 1 | Tris buffer 20 mmol/L, PEG 8000, pH 8.3. |
| Diluent | Preservative. |
| R 2 | Goat serum, anti-human α ₁ -antitrypsin, pH 7.5. |
| Antibody | Preservative. |
| Optional | Ref: 1102003 PROT CAL. |

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the α₁-antitrypsin calibrator by the corresponding factor stated in table below to obtain the α₁-antitrypsin concentration of each dilution.

| Calibrator dilution | 1 | 2 | 3 | 4 | 5 | 6 |
|---------------------|-----|-----|------|-----|------|-----|
| Calibrator (μL) | -- | 10 | 25 | 50 | 75 | 100 |
| NaCl 9 g/L (μL) | 100 | 90 | 75 | 50 | 25 | - |
| Factor | 0 | 0.1 | 0.25 | 0.5 | 0.75 | 1.0 |

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

REFERENCE VALUES

Newborn: Between 124 - 348 mg/dL.

Adults: 90 - 200 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINLAB 180 APPLICATION

| | | | |
|---------------------|-------------|------------------|-----------|
| Name | Antitrypsin | Ref. male low | 90 mg/dL |
| Abbr. Name | ATRYP | Ref. male high | 200 mg/dL |
| Mode | Twopoints | Ref. female low | 90 mg/dL |
| Wavelength | 340 nm | Ref. female high | 200 mg/dL |
| Units | mg/dL | Ref. Ped. Low | 124 mg/dL |
| Decimals | 0 | Ref. Ped. High | 348 mg/dL |
| Low Conc. | 16 mg/dL | Control 1 | * |
| High Conc. | 400 mg/dL | Control 2 | * |
| Calibrator name | SP CAL | Control 3 | * |
| Prozone check | No | Correlat. factor | 1.000 |
| | | Correlat. offset | 0.000 |
| DUAL MODE | | | |
| Sample blank | No | | |
| R1 bottle (mL) | 25 mL | | |
| normal volume | 318 μL | | |
| rerun volume | 318 μL | | |
| Sample | | | |
| normal volume | 2.0 μL | | |
| rerun volume | 2.0 μL | | |
| R2 bottle (mL) | 5 mL | | |
| normal volume | 80 μL | | |
| rerun volume | 80 μL | | |
| Predilución | No | | |
| Slope blank | No | | |
| Point one,two | -3.130 sec. | | |
| Reagent blank | No | | |
| Low Absorbance | -0.100 Abs | | |
| High Absorbance | 3.000 Abs | | |
| R. Abs. L. Limit | -0.100 Abs | | |
| R. Abs. H. Limit | 3.000 Abs | | |
| Substr.Depletion | 3.000 Abs | | |

PERFORMANCE CHARACTERISTICS

Measurement range: Up to 500 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

Detection Limit: Values less than 16 mg/dL give non-reproducible results.

Prozone effect: No prozone effect was detected upon 750 mg/dL

Precision: The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

| EP5 | CV (%) | | |
|-------------|-------------|------------|-------------|
| | 34.39 mg/dl | 92.7 mg/dl | 181.8 mg/dl |
| Total | 4.2% | 2.6% | 2.8% |
| Within Run | 0.8% | 1.1% | 1.6% |
| Between Run | 3.8% | 2.4% | 2.3% |
| Between Day | 1.6% | 0% | 0% |

Sensitivity: Δ 3.4 mA / mg/dL.

Accuracy: Results obtained using this reagent (y) were compared to those obtained using the Beckman Array 360 CE. 35 samples ranging from 70 to 250 mg/dL of α₁antitrypsin were assayed. The correlation coefficient (r) was 0.92 and the regression equation y = 0.8567x +26.5.

The results of the performance characteristics depend on the used analyzer.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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| EP5 | CV (%) | | |
|--------------|-------------|------------|-------------|
| | 34.39 mg/dl | 92.7 mg/dl | 181.8 mg/dl |
| Total | 4.2% | 2.6% | 2.8% |
| Intra-series | 0.8% | 1.1% | 1.6% |
| Entre-series | 3.8% | 2.4% | 2.3% |
| Entre-días | 1.6% | 0% | 0% |

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BIBLIOGRAFIA

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3tn ed. AACC Pres, 1997.

PRESENTACION

Ref: SP1102054

Cont.

R1: 2 x 20 mL

R2: 2 x 5 mL