# α<sub>1</sub>-Acid Glycoprotein

Turbidimetry

## Quantitative determination of α<sub>1</sub>-Acid Glycoprotein (α<sub>1</sub>-Ac GLY) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Anti- α<sub>1</sub>-Ac Gly antibodies when mixed with samples containing α<sub>1</sub>-Ac Gly, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α<sub>1</sub>-Ac Gly concentration of the patient sample, that can be quantified by comparison from a calibrator of known α<sub>1</sub>-Ac Gly concentration.

### CLINICAL SIGNIFICANCE

α<sub>1</sub>-Ac Glycoprotein (also known as orosomucoid) is a glycoprotein synthesized by hepatic parenchymal cells, but granulocytes and monocytes may also contribute significantly to plasma levels in sepsis. It has long been known to bind a large number of basic and lipophilic compounds (progesterone and related hormones).

It is an acute phase response protein that shows a 3 to 4-fold increase in most conditions associated with inflammation or tissue necrosis, and may be one of the most reliable indicators of clinical activity of ulcerative colitis. Levels also are increased by glucocorticoid effect.

Synthesis and plasma levels are decreased by estrogens.

### REAGENTS

<b>R 1</b> Diluent	Tris buffer 20 mmol/L, PEG 8000, pH, 8,3. Sodium azide 0,95 g /L.
<b>R 2</b> Antibody	Goat serum, anti-human α <sub>1</sub> -Ac Glycoprotein, pH 7.5. Sodium azide 0.95 g/L.
<b>Optional</b>	Ref: 1102003 PROT CAL.

### CALIBRATION

The assay has been standardized against the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM). It must be used the PROT CAL to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

### PREPARATION

**Reagents:** Ready to use.

**Calibration Curve:** Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the α<sub>1</sub>-Ac Glycoprotein calibrator by the corresponding factor stated in table below to obtain the α<sub>1</sub>-Ac Glycoprotein concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (μL)	--	10	25	50	75	100
NaCl 9 g/L (μL)	100	90	75	50	25	-
Factor	0	0.1	0,25	0,5	0,75	1,0

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

### ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

### REFERENCE VALUES

Between 50 - 120 mg/dL. Each laboratory should establish its own reference range.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## SPINLAB 180 APPLICATION

Name	Glycoprot.	Ref. male low	50 mg/dL
Abbr. Name	GLY	Ref. male high	120 mg/dL
Mode	Twopoints	Ref. female low	50 mg/dL
Wavelength	340 nm	Ref. female high	120 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	10 mg/dL	Control 1	*
High Conc.	250 mg/dL	Control 2	*
Calibrator name	SP CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
<b>Sample blank</b>	<b>No</b>		
R1 bottle (mL)	25 mL		
normal volume	240 μL		
rerun volume	240 μL		
Sample			
normal volume	2.0 μL		
rerun volume	2.0 μL		
R2 bottle (mL)	5 mL		
normal volume	60 μL		
rerun volume	60 μL		
Predilución	No		
Slope blank	No		
Point one,two	6.236 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

### PERFORMANCE CHARACTERISTICS

**Measurement range:** Up to 250 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

**Detection Limit:** Values less than 12.9 mg/dL give non-reproducible results.

**Prozone effect:** No prozone effect was detected upon 1000 mg/dL.

**Sensitivity:** Δ 5.0 mA / mg/dL.

**Precision:** The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	56.4 mg/dl	112.07 mg/dl
Total	3.3%	3.1%
Within Run	1.1%	1.6%
Between Run	3%	2.1%
Between Day	0.7%	1.6%

**Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the method Immage from Beckman. 51 samples ranging from 50 to 120 mg/dL of α<sub>1</sub>-Ac Gly were assayed. The correlation coefficient (r) was 0.95 and the regression equation y = 0.9304x + 6.5367.

The results of the performance characteristics depend on the used analyzer.

### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Bienvenue J et al. Clin Chem Clin Biochem 1981; 27: 721-726.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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