**Quantitative determination of Transferrin (TRF)**
**IVD**

Store 2 - 8°C.

**PRINCIPLE OF THE METHOD**

Anti-transferrin antibodies when mixed with samples containing TRF, form insoluble complexes. These complexes cause an absorbance change, dependent upon the TRF concentration of the patient sample, that can be quantified by comparison from a calibrator of known TRF concentration.

**CLINICAL SIGNIFICANCE**

Transferrin is a plasma protein that contains a single polypeptide chain with approximately 6% carbohydrate. It is synthesized in the liver and transfers iron through the serum.

Evaluation of plasma TRF levels is useful for the differential diagnosis of anemia and for monitoring its treatment. In the hypochromic anemia (iron deficiency), the TRF level is increased. On the other hand, if the anemia is due to a failure to incorporate iron into erythrocytes, the TRF level is normal or low but the protein is highly saturated with iron. In iron overload, the TRF concentration is normal but saturation exceeds 55% and may be as great as 90%.

TRF concentration may, in fact, be used for assessing nutritional status. In the congenital defect atransferrinemia, very low level of TRF is accompanied by iron overload and a severe hypochromic anemia that is resistant to iron therapy. High levels of TRF occur in pregnancy and during estrogen administration.

**REAGENTS**

|                 |                                             |
|-----------------|---------------------------------------------|
| <b>R 1</b>      | Tris buffer 20 mmol/L, PEG 8000, pH 8.3.    |
| Diluent         | Preservative.                               |
| <b>R 2</b>      | Goat serum, anti-human transferrin, pH 7.5. |
| Antibody        | Preservative.                               |
| <b>Optional</b> | Ref: 1102003 PROT CAL.                      |

**CALIBRATION**

The assay has been standardized against the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM). It must be used the PROT CAL to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

**PREPARATION**

**Reagents:** Ready to use.

**Calibration Curve:** Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the transferrin calibrator by the corresponding factor stated in table below to obtain the transferrin concentration of each dilution:

| Calibrator dilution | 1   | 2   | 3    | 4   | 5    | 6   |
|---------------------|-----|-----|------|-----|------|-----|
| Calibrator (µL)     | --  | 10  | 25   | 50  | 75   | 100 |
| NaCl 9 g/L (µL)     | 100 | 90  | 75   | 50  | 25   | -   |
| Factor              | 0   | 0.1 | 0.25 | 0.5 | 0.75 | 1.0 |

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

**ADDITIONAL EQUIPMENT**

- Spinlab 180autoanalyzer
- Laboratory equipment

**SAMPLES**

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**REFERENCE VALUES**

Between 200 - 360 mg/dL. Each laboratory should establish its own reference range.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**SPINLAB 180 APPLICATION**

|                     |             |                  |           |
|---------------------|-------------|------------------|-----------|
| Name                | Transferrin | Ref. male low    | 200 mg/dL |
| Abbr. Name          | TRF         | Ref. male high   | 360 mg/dL |
| Mode                | Twopoints   | Ref. female low  | 200 mg/dL |
| Wavelength          | 340 nm      | Ref. female high | 360 mg/dL |
| Units               | mg/dL       | Ref. Ped. Low    | *         |
| Decimals            | 0           | Ref. Ped. High   | *         |
| Low Conc.           | 1 mg/dL     | Control 1        | *         |
| High Conc.          | 750 mg/dL   | Control 2        | *         |
| Calibrator name     | SP CAL      | Control 3        | *         |
| Prozone check       | No          | Correlat. factor | 1.000     |
|                     |             | Correlat. offset | 0.000     |
| <b>DUAL MODE</b>    |             |                  |           |
| <b>Sample blank</b> | <b>No</b>   |                  |           |
| R1 bottle (mL)      | 25 mL       |                  |           |
| normal volume       | 240 µL      |                  |           |
| rerun volume        | 240 µL      |                  |           |
| Sample              |             |                  |           |
| normal volume       | 2.0 µL      |                  |           |
| rerun volume        | 2.0 µL      |                  |           |
| R2 bottle (mL)      | 5 mL        |                  |           |
| normal volume       | 60 µL       |                  |           |
| rerun volume        | 60 µL       |                  |           |
| Predilución         | No          |                  |           |
| Slope blank         | No          |                  |           |
| Point one,two       | -3.130 sec. |                  |           |
| Reagent blank       | No          |                  |           |
| Low Absorbance      | -0.100 Abs  |                  |           |
| High Absorbance     | 3.000 Abs   |                  |           |
| R. Abs. L. Limit    | -0.100 Abs  |                  |           |
| R. Abs. H. Limit    | 3.000 Abs   |                  |           |
| Substr.Depletion    | 3.000 Abs   |                  |           |

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**PERFORMANCE CHARACTERISTICS**

**Measurement range:** Up to 750 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

**Detection Limit:** 0 mg/dL.

**Prozone effect:** No prozone effect was detected upon 2000 mg/dL.

**Sensitivity:** Δ 3.0 mA / mg/dL (94 mg/dL).

**Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

| EP5         | CV (%) |             |              |
|-------------|--------|-------------|--------------|
|             |        | 77.02 mg/dl | 206.99 mg/dl |
| Total       | 5.4%   | 2.5%        | 5.4%         |
| Within Run  | 1%     | 0.8%        | 1.2%         |
| Between Run | 1.7%   | 1.3%        | 2.1%         |
| Between Day | 5%     | 2%          | 4.9%         |

**Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Immage from Beckman. 100 samples ranging from 50 to 700 mg/dL of TRF were assayed. The correlation coefficient (r) was 0.95 and the regression equation  $y = 1.046x + 3.843$ .

The results of the performance characteristics depend on the used analyzer.

**BIBLIOGRAPHY**

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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### BIBLIOGRAFIA

1. McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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|                |              |                   |           |
|----------------|--------------|-------------------|-----------|
| Nome           | Transferrina | Ref. Homem Inf.   | 200 mg/dL |
| Nome abreviado | TRF          | Ref. Homem Sup.   | 360 mg/dL |
| Modo           | Twopoint     | Ref. Mulher Inf.  | 200 mg/dL |
| Comp. ondas    | 340 nm       | Ref. Mulher Sup.  | 360 mg/dL |
| Unidades       | mg/dL        | Ref. Ped. Inf.    | *         |
| Decimais       | 0            | Ref. Ped. Sup.    | *         |
| Conc. Inferior | 1 mg/dL      | Control 1         | *         |
| Conc. Superior | 750 mg/dL    | Control 2         | *         |
| Calibrador     | CAL PS       | Control 3         | *         |
| Check prozona  | Não          | Factor correl.    | 1.000     |
|                |              | Offset de correl. | 0.000     |

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| EP5         | CV (%)      |              |           |
|-------------|-------------|--------------|-----------|
|             | 77,02 mg/dL | 206.99 mg/dL | 377 mg/dL |
| Total       | 5,4%        | 2,5%         | 5,4%      |
| Within Run  | 1%          | 0,8%         | 1,2%      |
| Between Run | 1,7%        | 1,3%         | 2,1%      |
| Between Day | 5%          | 2%           | 4,9%      |

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### BIBLIOGRAFIA

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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Ref.: SP1102134

Cont.

R1: Diluyente: 2 x 20 mL

R2: Anticorpo: 2 x 5 mL