Antithrombin III

Turbidimetry

Quantitative determination of Antithrombin III (ATHROM-III) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti-antithrombin III antibodies when mixed with samples containing antithrombin-III, form insoluble complexes. These complexes cause an absorbance change, dependent upon the antithrombin-III concentration of the patient sample, that can be quantified by comparison from a calibrator of known antithrombin-III concentration.

CLINICAL SIGNIFICANCE

Antithrombin III is a protein synthesized in the liver, normally present in the human plasma. It is the major inhibitor of the thrombin, and inhibits coagulation and limits the forming of blood clots. Antithrombin-III is also capable of activating other components of the coagulation cascade (eg, factor Xa), as well as plasmin.

Antithrombin-III deficiency can cause or lead to thrombosis, a clot forming in a blood vessel. Clots forming in the legs and pulmonary embolism are most commonly reported. Antithrombin-III deficiency is usually inherited and affects males and females equally. All family members should be tested if there is history of the disease.

Acquired antithrombin-III deficiency can occur as a result of other conditions. It has been reported in patients with liver diseases, patients receiving certain kinds of chemotherapy, and patients using oral contraceptives.

REAGENTS

R 1	Tris buffer 20 mmol/L, PEG 8000, pH, 8,3.
Diluent	Preservative.
R 2	Goat serum, anti-human antithrombin III, pH 7.5.
Antibody	Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

It must be used the PROT CAL to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every week, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the ATHROM-III calibrator by the corresponding factor stated in table below to obtain the antithrombin-III concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. Sodium citrate should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 17 - 30 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINLAB 180 APPLICATION

Name	Antithrombin III	Ref. male low	17.00 mg/dL
Abbr. Name	AT-III	Ref. male high	30.00 mg/dL
Mode	Twopoints	Ref. female low	17.00 mg/dL
Wavelength	340 nm	Ref. female high	30.00 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	7.40 mg/dL	Control 1	*
High Conc.	70.00 mg/dL	Control 2	*
Calibrator name	SP CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	6.0 µL		
rerun volume	6.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	-3.236 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

PERFORMANCE CHARACTERISTICS

1. Measurement range: Up to 70 mg/dL (Nota 1), under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Limit detection: Values less than 7.4 mg/dL give non-reproducible results.

3. Prozone effect: No prozone effect was detected upon 200 mg/dL

4. Sensitivity: Δ 7.5 mA / mg/dL.

5. Precision: The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	17.85 mg/dl	35.93 mg/dl
Total	3.2%	3.5%
Within Run	0.8%	1%
Between Run	2.4%	2.4%
Between Day	2%	2.3%

6. Accuracy: Results obtained using this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
3. Buller HR et al. Critical care Medicine 1982; 10: 311.
4. Kauffman et al. Am J Med 1978; 65: 607.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Buller HR et al. Critical care Medicine 1982; 10: 311.
- Kauffman et al. Am J Med 1978; 65: 607.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

PRESENTACION

Ref: SP1102014 Cont. R1: 2 x 20 mL
R2: 2 x 5 mL