Transferrin

Turbidimetry

Quantitative determination of Transferrin (TRF) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-transferrin antibodies when mixed with samples containing TRF, form insoluble complexes. These complexes cause an absorbance change, dependent upon the TRF concentration of the patient sample, that can be quantified by comparison from a calibrator of known TRF concentration.

CLINICAL SIGNIFICANCE

Transferrin is a plasma protein that contains a single polypeptide chain with approximately 6% carbohydrate. It is synthesized in the liver and transfers iron through the serum.

Evaluation of plasma TRF levels is useful for the differential diagnosis of anemia and for monitoring its treatment. In the hypochromic anemia (iron deficiency), the TRF level is increased. On the other hand, if the anemia is due to a failure to incorporate iron into erythrocytes, the TRF level is normal or low but the protein is highly saturated with iron. In iron overload, the TRF concentration is normal but saturation exceeds 55% and may be as great as 90%.

TRF concentration may, in fact, be used for assessing nutritional status. In the congenital defect atransferrinemia, very low level of TRF is accompanied by iron overload and a severe hypochromic anemia that is resistant to iron therapy. High levels of TRF occur in pregnancy and during estrogen administration.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Sodium azide 0.95 g/L.
R 2 Antibody	Goat serum, anti-human transferrin, pH 7.5. Sodium azide 0.95 g/L.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM). It must be used the PROT CAL to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL Calibrator dilutions in CINA 9 g/L as diluent. Multiply the concentration of the transferrin calibrator by the corresponding factor stated in table below to obtain the transferrin concentration of each dilution:

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SPINLAB 180 APPLICATION

Name	Transferrin	Ref. male low	200 mg/dL
Abbr. Name	TRF	Ref. male high	360 mg/dL
Mode	Twopoints	Ref. female low	200 mg/dL
Wavelength	340 nm	Ref. female high	360 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	1 mg/dL	Control 1	*
High Conc.	750 mg/dL	Control 2	*
Calibrator name	SP CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	2.0 µL		
rerun volume	2.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	-3.130 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 200 - 360 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. *McNeely Clinical Guide to Laboratory Tests*, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. *Dati F et al. Eur J Clin Chem Clin Biochem* 1996; 34:517-520.
3. *Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.*
4. *Kreutzer HJH. J Clin Chem Clin Biochem* 1976; 14: 401-406
5. *Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.*
6. *Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.*

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