Lp(a)-turbilatex

Latex turbidimetry

Quantitative determination of Lipoprotein (a) (Lp(a)) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The Lp(a)-turbilatex is a quantitative turbidimetric test for the measurement of Lp(a) in human serum or plasma. Latex particles coated with antibodies anti-Lp(a) are agglutinated when mixed with samples containing Lp(a). The agglutination causes an absorbance change, dependent upon the Lp(a) contents of sample that can be quantified by comparison from a calibrator of known Lp(a) concentration.

CLINICAL SIGNIFICANCE

Lp(a) is a low density lipoprotein-like particle containing apolipoprotein B-100 disulphide-linked to one large glycoprotein called apolipoprotein (a). Many investigators have confirmed that a high Lp(a) concentration represents an indicator of risk for cardiovascular disease, especially when serum LDL-cholesterol or Apo B are elevated. The quantification of Lp(a) in serum or plasma is important for identification of individuals at risk for developing atherosclerosis.

REAGENTS

| | |
|---------------------|---|
| Diluent (R1) | Glycine buffer 50 mmol/L, pH 9.0. Sodium azide 0.95 g/L. |
| Latex (R2) | Latex particles coated with mouse monoclonal anti-human Lp(a), pH 8.2. Sodium azide 0.95 g/L. |
| Optional | Ref.: 1107022 Lp(a) Calibrator. Ref.: 1107024 Lp(a) Control. |

CALIBRATION

The sensitivity of the assay and the target value of the calibrator have been standardized against an Internal Reference Material. It is not recommended the use of other commercially available Lp(a) calibrators.

PREPARATION

Calibration Curve Prepare the following Lp(a) calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the Lp(a) calibrator by the corresponding factor stated in table below to obtain the Lp(a) concentration of each dilution.

| Calibrator dilution | 1 | 2 | 3 | 4 | 5 |
|-----------------------|-----|-------|------|-----|-----|
| Lp(a) Calibrator (µL) | -- | 12.5 | 25 | 50 | 100 |
| NaCl 9 g/L (µL) | 100 | 87.5 | 75 | 50 | - |
| Factor | 0 | 0.125 | 0.25 | 0.5 | 1.0 |

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: Presence of particles and turbidity.

Do not freeze; frozen latex and diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Normal values up to 30 mg/dL. Each laboratory should establish its own reference range.

SPINLAB 180 APPLICATION

| | | | |
|-----------------|------------|------------------|------------|
| Name | Lp(a) | Ref. male low | 0.0 mg/dL |
| Abbr. Name | Lp(a) | Ref. male high | 30.0 mg/dL |
| Mode | Twopoints | Ref. female low | 0.0 mg/dL |
| Wavelength | 578 nm | Ref. female high | 30.0 mg/dL |
| Units | mg/dL | Ref. Ped. Low | * |
| Decimals | 1 | Ref. Ped. High | * |
| Low Conc. | 12.0 mg/dL | Control 1 | * |
| High Conc. | 80.0 mg/dL | Control 2 | * |
| Calibrator name | CAL Lp(a) | Control 3 | * |
| Prozone check | No | Correlat. factor | 1.000 |
| | | Correlat. offset | 0.000 |

DUAL MODE

| | |
|----------------|-------------|
| Sample blank | No |
| R1 bottle (mL) | 25 mL |
| normal volume | 240 µL |
| rerun volume | 240 µL |
| Sample | |
| normal volume | 5.0 µL |
| rerun volume | 3.0 µL |
| R2 bottle (mL) | 5 mL |
| normal volume | 60 µL |
| rerun volume | 60 µL |
| Predilución | No |
| Slope blank | No |
| Point one,two | 6, 236 sec. |

| | |
|-------------------|------------|
| Reagent blank | Yes |
| Low Absorbance | -0.100 Abs |
| High Absorbance | 3.000 Abs |
| R. Abs. L. Limit | -0.100 Abs |
| R. Abs. H. Limit | 3.000 Abs |
| Substr. Depletion | 3.000 Abs |

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Lp(a) Control Ref.: 1107024.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

Hemoglobin (5 g/L), bilirubin (20 mg/dL), plasminogen (680 mg/dL), ascorbic ac. (200 mg/dL), rheumatoid factors (100 IU/mL) and lipemia (20 g/L), do not interfere. Other substances may interfere⁵.

NOTES

1. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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