Cholinesterase bi



Butyrylthiocholine. Kinetic

Quantitative determination of cholinesterase (CHE) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Cholinesterase hydrolysed butyrylthiocholine to butyrate and thiocholine. Thiocholine reacts with 5,5'-dithiobis-2-nitrobenzoic acid (DTNB) to form nitro-2-mercapto-5-benzoato according the following reactions:

The rate of nitro-2-mercapto-5-benzoato formation, measured photometrically, is proportional to the enzymatic activity of cholinesterase in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Cholinesterase is an enzyme present in plasma and synthesized by the liver. Its true physiological function is unknown, so its function may be to hydrolyze choline in plasma. Cholinesterase activity is usually measured for liver function, is a sensitive test of exposure to pesticides organophosphorus and identification of patients with the atypical form of enzyme whose presents high sensitivity to succinylcholine^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Phosphate pH 7.7	50 mmol/L
Buffer	5,5 DTNB	0.25 mmol/L
R 2	Butyrylthiocholine	7 mmol/L
Substrate		

PREPARATION

- Working reagent (WR₁):
Dissolve (→) the contents of each vial of R 1 in 50 mL of distilled water.
 - Working reagent (WR₂):
Dissolve (→) the contents of each vial of R 2 in 3 mL of distilled water.
- Cap vial and mix gently to dissolve contents.
Stability: 8 weeks at 2-8°C. Avoid direct sunlight.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.
Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm ≥ 1.20.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 405 nm.
- Thermostatic bath at 25°C, 30°C o 37°C (± 0.1°C)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or heparinized or EDTA plasma¹: Stability 7 days at 2-8°C.

PROCEDURE

1. Assay conditions:
Wavelength: 405 nm
Cuvette: 1 cm light path
Constant temperature 25°C/30°C / 37°C
2. Adjust the instrument to zero with distilled water or air.
3. Pipette into a cuvette:

	25 - 30°C	37°C
WR ₁ (mL)	1.5	1.5
WR ₂ (μL)	50	50
Sample (μL)	10	--
Sample diluted 1/2 with NaCl 9 g/L (μL)	--	10

4. Mix and wait 30 seconds.

5. Read initial absorbance (A) of the sample, start the stopwatch and read the absorbance at 30 seconds intervals thereafter for 1.5 minutes.
6. Calculate the difference between absorbances and the average absorbance differences per 30 seconds (ΔA/30 s).

CALCULATIONS

$$25^\circ - 30^\circ\text{C} \quad \Delta A / 30 \text{ s} \times 23460 = \text{U/L}$$

$$37^\circ\text{C} \quad \Delta A / 30 \text{ s} \times 46920 = \text{U/L}$$

Calculating factor in automatic analyzers (ΔA/min.) is 23460 at 37°C.

Units: One international unit (IU) is the amount of enzyme that transforms 1 μmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.24	1.55
30°C	0.81	1.00	1.26
37°C	0.64	0.80	1.00

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

25°C 30°C 37°C

3000 – 9300 U/L 3714 – 11513 U/L 4659 – 14443 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 6 U/L to *linearity limit* of 9000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/5 with NaCl 9 g/L and multiply the result by 5.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (U/L)	6571	3309	6557	3342
SD	45.7	24.6	53.4	43.5
CV (%)	0.70	0.74	0.81	1.30

Sensitivity: 1 U/L = 0.00006 ΔA / 30 s.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Moderate hemolysis will not interfere in the results.^{1,2}. A list of drugs and other interfering substances with cholinesterase determination has been reported by Young et al.^{3,4}.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

1. King M. Cholinesterase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1108-1111.
2. Whittaker M. et al. Comparison of a Commercially Available Assay System with Two Reference Methods for the Determination of Plasma Cholinesterase Variants. Clin. Chem 1983;(29/10); 1746-1760.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAPHIE

- King M. Cholinesterase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1108-1111.
- Whittaker M. et al. Comparation of a Commercially Available Assay System with Two Reference Methods for the Determination of Plasma Cholinesterase Variants. Clin. Chem 1983; (29/10): 1746-1760.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

- King M. Cholinesterase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1108-1111.
- Whittaker M. et al. Comparasion of a Commercially Available Assay System with Two Reference Methods for the Determination of Plasma Cholinesterase Variants..Clin. Chem 1983;(29/10); 1746-1760.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: 1001102	Cont.	R1: 4 → 50 mL
		R2: 4 → 3 mL