# Glucose-6-phosphate dehydrogenase

Enzymatic-UV

## Quantitative determination of G-6-PDH IVD


Store at 2-8°C

### INTENDED USE

For the quantitative *in vitro* determination of Glucose-6-Phosphate Dehydrogenase in erythrocytes.

### PRINCIPLE OF THE METHOD<sup>1,2,3</sup>

The enzyme activity is determined by measurement of the rate of absorbance change at 340 nm due to the reduction of NADP<sup>+</sup>.

### CLINICAL SIGNIFICANCE

G-6PDH deficiency is one of the most common human enzyme deficiencies, and it is estimated to affect more than 400 million people worldwide. Although most of the enzyme-deficient subjects are asymptomatic, deficient individuals may show episodic hemolytic anemia induced by infections or certain drugs and a spontaneous chronic nonspherocytic hemolytic anemia.

### REAGENTS

<b>R 1</b> Buffer	Triethanolamine pH 7,6 EDTA	31,7 mmol/L 3,2 mmol/L
<b>R 2</b>	NADP	0,34 mmol/L
<b>R3</b> Substrate	Glucose-6-phosphate	0,58 mmol/L
<b>R4</b>	Digitonin	
<b>OPTIONAL</b>	G-6-PDH 2 levels control ref.1002520	

### PRECAUTIONS

- **R2, R3, Control:** H302-Harmful if swallowed. H412-Harmful to aquatic life with long lasting effects.

Follow the precautionary statements given in MSDS and label of the product.

### PREPARATION

- **R1 - R4** Ready to use.

- **R2 - R3** Reconstitute the contents of each bottle with 2 mL of redistilled water.

### STORAGE AND STABILITY

R1 and R4 are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. R2 and R3 are stable for 4 weeks at 2-8°C. Do not use reagents over the expiration date.

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 37° C (± 0,1°C)
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment <sup>(Note 1)</sup>.

### SAMPLES

- Erythrocytes: Preparation: Wash 0,2 ml of blood with 2 mL aliquots of 0,9% NaCl solution. Centrifuge after each wash for 10 min at around 3000 rpm. Repeat 3 times. Suspend the washed centrifuged erythrocytes in 0,5 mL of R4 and let stand for 15 min at +4°C and then centrifuge again. Use the supernatant in the assay within 2 hours.

### PROCEDURE

- Assay conditions:  
Wavelength: .....340nm  
Cuvette: .....1 cm light path  
Constant temperature .....37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a test tube <sup>(Note 2)</sup>.

	Sample
Reagent 1 (µL)	1000
Reagent 2 (µL)	30
Haemolysate (µL)	15

4. Mix, incubate at 37°C for 5 minute.

5. Add:

	Sample
Reagent 3 (µL)	15

6. Mix, read initial absorbance and start timer simultaneously. Read again after 1, 2 and 3 minutes.

### CALCULATIONS

To calculate the G-6-PDH activity use the following calculations:  
 mU/erythrocytes per mL blood\* = 33650 x ΔA 340 nm/min  
 mU/erythrocytes per mL blood\* = 60571 x Δ A Hg 365 nm/min  
 mU/erythrocytes per mL blood\* = 34304 x Δ A Hg 334 nm/min

\*G-6-PDH activity is expressed as mU/10<sup>9</sup> erythrocytes or as mU/g haemoglobin. To calculate G-6-PDH activity as mU/10<sup>9</sup> erythrocytes for comparison with the normal value, divide the calculated activity (mU/erythrocytes per mL blood) with the RBC's count per mL.

To calculate G-6-PDH activity as mU/g haemoglobin use the following equation:

$$\text{G-6-PDH mU/gHb} = \frac{\text{mU.erythrocytes per mL} \times 100}{\text{Hb (g/dL)}}$$

$$\begin{aligned} 100 &= \text{Factor to convert from mL to dL} \\ \text{Hb (g/dL)} &= \text{Haemoglobin concentration determined for each specimen} \end{aligned}$$

### Temperature Correction

Please note that the below temperature correction may be used for patient samples only.

When the temperature is 37°C, no temperature correction factor (TCF) is required. If results for patient samples are reported at a temperature other than 37°C, a TCF must be used.

Cuvette Temperature (°C)	TCF
25	2.076
30	1.515

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: G-6-PDH 2 levels control ref.1002520.

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

In erythrocytes: 245 - 299 mU/10<sup>9</sup> erythrocytes (37°C).

6,97 – 20,5	U/g Hb (37°C)	Blood Haemoglobin (g/dL)
		Adult Males 13 - 18
		Adult Females 11 - 16
		Newborns 14 - 23

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 154 U/L to *linearity limit* of 4000 U/L. If the results obtained were greater than linearity limit, using 0,2 mL of haemolysate with 1,8 mL of NaCl 9g/L and repeat the assay. Then multiply the result by 10.

#### Precision:

	Intraserie (n=20)		Interserie (n=20)	
Media (U/L)	730	1460	730	1460
SD	30,95	70,52	47,09	76,61
CV (%)	4,24	4,83	6,45	5,11

**Sensitivity:** 1 U/L = 0,000077 (A)

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,9807.

Regression equation: y = 1,0069x + 47,644.

The results of the performance characteristics depend on the analyzer used.

### INTERFERENCES

Reticulocytes have higher G-6-PDH levels than mature red cells. Therefore, it is not recommended that this assay be performed after a severe hemolytic crisis, since G-6-PDH levels may appear falsely elevated. Testing may be more helpful after the level of mature red cells have returned to normal.

### NOTES

- In order to avoid contamination it is recommended to use disposable material.
- Use clean disposable pipette for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers.**

### BIBLIOGRAPHY

- Kornberg, A. *et al.*, Methods in Enzymology 1, Academic Press, New York, 1955; p.323.
- Makarem, A., Clinical Chemistry-Principles and Techniques. 2<sup>nd</sup> Ed. R.F. Henry, D. C. Cannon, J.W. Winkelman, Editors. Harper and Row, Hagerstown [MD], 1974; 1128-1135.
- Lohr GW, Waller HD: Glucose-6-Phosphate Dehydrogenase. Methods of Enzymatic Analysis, 3<sup>rd</sup> Edition - Verlag Chemie, Wehnhheim: 1974; p. 636.

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