Cholinesterase

Butyrylthiocholine. Kinetic

Quantitative determination of cholinesterase (CHE) in serum and plasma IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The method uses butyrylthiocholine as the specific substrate for cholinesterase (CHE). Cholinesterase catalyses the hydrolysis of butyrylthiocholine substrate forming butyrate and thiocholine. Thiocholine reduces hexacyanoferrate (III) to hexacyanoferrate (II). The decrease in absorbance is directly proportional to CHE activity in the sample^{7,8,9}.

CLINICAL SIGNIFICANCE

Cholinesterase is an enzyme present in plasma and synthesized by the liver. Its true physiological function is unknown, so its function may be to hydrolyze choline in plasma. Cholinesterase activity is usually measured for liver function, is a sensitive test of exposure to pesticides organophosphorus and identification of patients with the atypical form of enzyme whose presents high sensitivity to succinyl-choline^{1,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Pyrophosphate pH 7,6	92 mmol/L
Buffer	Hexacyanoferrate (III)	2,5 mmol/L
R 2	Butyrylthiocholine	91 mmol/L
Substrate		

PREPARATION

Reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Stability: 30 days at 2-8°C after opening, if contamination avoided.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 405 nm.
- Thermostatic bath at 37°C (± 0,1°C)
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES^{1,2,3,4,5}

Fresh serum, plasma (EDTA heparin) not haemolyzed and promptly separated from the red blood cells. Do not use sodium fluoride as an anticoagulant because it inhibits cholinesterase. Stability: 7 days at 2-8°C or 3 months at -20°C.

PROCEDURE

- Assay conditions:
Wavelength: 405 nm
Cuvette: 1 cm light path
Temperature: 37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette ^(Note 1):

	Reagent blank	Calibrator / Sample
Reagent 1 (µL)	1000	1000
Distilled water (µL)	20	--
Calibrator / Sample (µL)	--	20

4. Mix and incubate for 5 minutes.

5. Add:

Reagent 2 (µL)	200	200
----------------	-----	-----

6. Mix carefully and read the absorbance A_1 after 90 seconds and read the second absorbance A_2 after other 90 seconds.

7. Calculate: $\Delta A = A_1 - A_2$.

CALCULATIONS

$$\text{CHE Activity} = \frac{\Delta A \text{ Sample} - \Delta A \text{ Blank}}{\Delta A \text{ Calibrator} - \Delta A \text{ Blank}} \times \text{Calibrator conc.} = \text{U/L}$$

Units: One international unit (IU) is the amount of enzyme that transforms 1 µmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

Conversion factor:

$$\text{CHE (U/L)} \times 0,001 = \text{CHE (KU/L)}$$

$$\text{CHE (U/L)} \times 0,01667 = \text{CHE (ukat/L)}$$

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Adults (37°C) Male: 5100-11700 U/L

Female: 4000-12600 U/L

In infants up to 6 months of age, cholinesterase activity is 40% to 50% higher than in adults. In young women, the enzyme activity is approximately 64% to 74% of that in adult males. The activity decreases during pregnancy.

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 150 U/L to linearity limit of 22000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/L)	SD	CV (%)	
Mean (mg/L)	7122	5182	6966	5060
SD	127,1	67,9	168,9	110,6
CV (%)	1,78	1,31	2,42	2,19

Sensitivity: 1 U/L = -0,000013 (A)

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)²: 0,998.

Regression equation: $y = 1,11 x + 384,58$.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

The test is not affected by the presence of ascorbic acid up to 20 mg/dL, haemoglobin up to 1000 mg/dL, bilirubin up to 60 mg/dL and lipids up to 1000 mg/dL.

NOTES

- In order to avoid contamination it is recommended to use disposable material.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- NCCLS Document: "Procedures for the Collection of Arterial Blood Specimens; Approved Standard - Third Edition (1999)".
- Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).
- EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of Good Laboratory Practice as specified in Council Directive 87/18/EEC.
- Jakobs, D.S., Kasten, B.L., Demmott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
- Deutsche Gesellschaft für Klinische Chemie. Proposal of Standard Methods for the determination of enzyme catalytic concentrations in serum and plasma at 37°C. II Cholinesterase (acylcholine acylhydrolase, EC 3.1.1.8). Eur.J.Clin.Chem.Clin. Biochem 30,163 (1992).

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- NCCLS Document: "Procedures for the Collection of Arterial Blood Specimens; Approved Standard - Third Edition (1999)".
- Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).
- EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of Good Laboratory Practice as specified in Council Directive 87/18/EEC.
- Jakobs, D.S., Kasten, Jr., B.L., Demmott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
- Deutsche Gesellschaft für Klinische Chemie. Proposal of Standard Methods for the determination of enzyme catalytic concentrations in serum and plasma at 37°C. II Cholinesterase (acetylcholine acylhydrolase, EC 3.1.1.8). Eur.J.Clin.Chem.Clin. Biochem 30,163 (1992).

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BIBLIOGRAPHIE

- NCCLS Document: "Procedures for the Collection of Arterial Blood Specimens; Approved Standard - Third Edition (1999)".
- Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).
- EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of Good Laboratory Practice as specified in Council Directive 87/18/EEC.
- Jakobs, D.S., Kasten, Jr., B.L., Demmott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
- Deutsche Gesellschaft für Klinische Chemie. Proposal of Standard Methods for the determination of enzyme catalytic concentrations in serum and plasma at 37°C. II Cholinesterase (acylcholine acylhydrolase, EC 3.1.1.8). Eur.J.Clin.Chem.Clin. Biochem 30,163 (1992).

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Réf. : 41210

Cont.

R1 : 1x100 mL R2 : 1x20 mL