Angiotensin Converting Enzyme

FAP Method

Quantitative determination of Angiotensin Converting Enzyme (ACE) in serum or plasma IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The ACE reagent is based on the FAP method. In this method the direct substrate N-[3-(2-furyl)-acryloyl]-L-phenylalanyl-glycylglycine (FAPGG) is hydrolysed to FAP and Glycylglycine. The hydrolysis of FAPGG by ACE results in a decrease in absorbance at 340nm.

CLINICAL SIGNIFICANCE

ACE reagent is for use in the determination of angiotensin converting enzyme (ACE) activity in serum or plasma at 340 nm. Angiotensin converting enzyme (ACE, dipeptidyl carboxypeptidase) is a glycoprotein peptidyl dipeptide hydrolase that cleaves histidylleucine dipeptide from angiotensin I, a relatively inactive decapeptide. The latter is converted to the potent vasoconstrictor, angiotensin II. ACE also inactivates bradykinin. Elevated levels of ACE activity occur in serum of patients with active sarcoidosis, and occasionally in premature infants with respiratory distress syndrome, in adults with tuberculosis, Gaucher's disease, leprosy, and in many other pathologic conditions involving lung and liver diseases^{6,7}.

REAGENTS

R	Goods buffer, pH 8.2	80 mmol/L
	FAPGG	0,50 mmol/L
	Sodium azide	< 0,1%

PREPARATION

The reagent is ready to use.

CALIBRATION

It is recommended to use the ACE Calibrator (Ref. 1002225).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C. After opening, the reagent is stable for 30 days when properly capped immediately after each opening and stored at 2-8°C protected from light.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 37° C (± 0.1°C)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Fresh serum or plasma (lithium or sodium-heparin) promptly separated from the red blood cells¹. Do not use: lipemic samples, haemolyzed samples or EDTA as an anticoagulant as it inhibits the ACE activity^{2,3,4,5,8}. Stability: 7 days at 2-8°C or 1 year at -20°C.

PROCEDURE

1. Assay conditions:
 Wavelength: 340nm(2nd 660nm)
 Cuvette: 1 cm light path
 Constant temperature: 37°C
 2. Adjust the spectrophotometer to zero with distilled water or air.
 3. Pipette into a cuvette:
- | | | |
|-------------|--------|------------|
| | Sample | Calibrator |
| R (mL) | 1,0 | 1,0 |
| Sample (µL) | 100 | |
| Calibrator | | 100 |
4. Mix and incubate 5 minutes at 37°C.
 5. Read initial absorbance of the sample and calibrator (A₁), start the stopwatch and read absorbance exactly after 5 minutes (A₂). Calculate the difference of absorbance for the sample and calibrator (A₂ - A₁).

CALCULATIONS

$$\frac{(A_2 - A_1) \text{ Sample}}{(A_2 - A_1) \text{ Calibrator}} \times \text{Conc. Calibrator} = \text{U/L of ACE in the sample}$$

Conversion factor

$$\text{ACE (U/L)} \times 0,001 = \text{ACE (kU/L)}$$

$$\text{ACE (U/L)} \times 0,01667 = \text{ACE (µkat/L)}$$

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: ACE Control (Ref.1002227). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹⁰

13,3 – 63,9 U/L
 These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: 5,4 - 160 U/L. If the sample concentration exceeds this value, dilute the sample 1:10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intraserie (n=20)		Interserie (n=20)	
Media (U/L)	3,71	111,9	40,1	114,7
SD	0,804	1,044	1,230	1,949
CV (%)	2,17	0,93	3,07	1,70

Sensitivity: 1U/L = 0,00041 (ΔA)

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results were the following:
 Correlation coefficient (r)²: 0,976.
 Regression equation: y = 0,98x - 0,56.
 The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Unconjugated Bilirubin up to 13 mg/dL, conjugated Bilirubin up to 26 mg/dL, haemoglobin up to 100 mg/dL and lipids up to 200 mg/dL. Captotril, and ACE inhibitory drug, used for the treatment of hypertension and some types congestive heart failure, will inhibit ACE activity in serum or plasma^{9,11}

NOTES

1.SPINREACT has instruction sheets for several automatic analyzers.

BIBLIOGRAPHY

1. NCCLS Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition (H3-A5). Wayne, PA: The National Clinical Laboratory Standards, 2003.
2. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne Pathogens.
3. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. Washington, DC: US Government Printing Office, January 2007.
4. World Health Organization. Laboratory Biosafety Manual, 3rd Ed. Geneva: World Health Organization, 2004.
5. Sewell DL, Bove KE, Callihan DR, et al. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Third Edition (M29-A3). Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
6. Pesce, A.J., Kaplan, L.A.: "Methods in Clinical Chemistry", Mosby Ed. (1987).
7. Burtis C.A., Ashwood E.R.: "Tietz Textbook of Clinical Chemistry", W.B. Saunders Company Ed. (3rd edition, 1999).
8. Guder W.G.: "The Quality of Diagnostic Sample". Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine. (1st Edition - 2001).
9. Jakobs, D.S., Kasten, Jr., B.L., DeMott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
10. Silvia Camós, M. Jesús Cruz, Ferran Morell and Esther Solé: "Genetic-based reference values for angiotensin converting enzyme (ACE) according to I/D polymorphism in a Spanish population sample". Clin Chem Lab Med 2012;50(10):1749 -1753.
11. Beneteau B, Baudin B, Morgant G, Giboudeau J, Baumann FC. Automated kinetic assay of angiotensin-converting enzyme in serum. Clin Chem 1986; 32: 884-6.

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BIBLIOGRAFIA

- NCCLS Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition (H3-A5). Wayne, PA: The National Clinical Laboratory Standards, 2003.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne Pathogens.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. Washington, DC: US Government Printing Office, January 2007.
- World Health Organization. Laboratory Biosafety Manual, 3rd Ed. Geneva: World Health Organization, 2004.
- Sewell DL, Bove KE, Callihan DR, et al. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Third Edition (M29-A3). Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
- Pesce, A.J., Kaplan, L.A.: "Methods in Clinical Chemistry", Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: "Tietz Textbook of Clinical Chemistry", W.B. Saunders Company Ed. (3rd edition, 1999).
- Guder W.G.: "The Quality of Diagnostic Sample". Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine. (1st Edition - 2001).
- Jakobs, D.S., Kasten, Jr., B.L., DeMott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
- Silvia Camós, M. Jesús Cruz, Ferran Morell and Esther Solé: "Genetic-based reference values for angiotensin-converting enzyme (ACE) according to I/D polymorphism in a Spanish population sample". Clin Chem Lab Med 2012;50(10):1749 -1753.
- Beneteau B, Baudin B, Morgant G, Giboudeau J, Baumann FC. Automated kinetic assay of angiotensin-converting enzyme in serum. Clin Chem 1986; 32: 884-6.

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Ref: 41205

Cont.

R: 2 x 50 mL