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
- Naito H.K. Cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.
- Meiattini F. et al. The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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


| | | |
|--------------|-------|---|
| Ref: 1001090 | | R1: 10 x 50 mL, R2: 10 → 50 mL, CAL: 1 x 5 mL |
| Ref: 1001091 | Cont. | R1: 10 x 20 mL, R2: 10 → 20 mL, CAL: 1 x 5 mL |
| Ref: 1001092 | | R1: 4 x 125 mL, R2: 4 → 125 mL, CAL: 1 x 5 mL |
| Ref: 1001093 | | R1: 4 x 250 mL, R2: 4 → 250 mL, CAL: 1 x 5 mL |

Quantitative determination of cholesterol IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The cholesterol present in the sample originates a coloured complex, according to the following reaction:

The intensity of the color formed is proportional to the cholesterol concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Cholesterol is a fat-like substance that is found in all body cells. The liver makes all of the cholesterol the body needs to form cell membranes and to make certain hormones.

The determination of serum cholesterol is one of the important tools in the diagnosis and classification of lipemia. High blood cholesterol is one of the major risk factors for heart disease^{5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| | | |
|---|--|------------|
| R 1 Buffer | PIPES pH 6,9 | 90 mmol/L |
| | Phenol | 26 mmol/L |
| R 2 Enzymes ^(Note 2) | Cholesterol esterase (CHE) | 300 U/L |
| | Cholesterol oxidase (CHOD) | 300 U/L |
| | Peroxidase (POD) | 1250 U/L |
| | 4 - Aminophenazone (4-AP) | 0,4 mmol/L |
| CHOLESTEROL CAL | Cholesterol aqueous primary standard 200 mg/dL Contains Triton X-114 10-15% | |

PRECAUTIONS

CAL: H225- Highly flammable liquid and vapour. H318- Causes serious eye damage. H412- Harmful to aquatic life with long lasting effects. Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Working reagent (WR): Dissolve (→) the contents of one vial R 2 Enzymes in one bottle of R 1 Buffer.

Cap and mix gently to dissolve contents.

(WR) is stable: 4 months at 2-8°C or 40 days at 15-25°C.

Avoid direct sunlight.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm ≥ 0,1.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 505 nm (500-550).
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma^{1,2}: Stability of the sample for 7 days at 2-8°C or freezing at -20°C will keep samples stable for a 3 months.

PROCEDURE

- Assay conditions:
Wavelength:..... 505 nm (500-550)
Cuvette:1 cm light path
Temperature:.....37°C /15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette^(Note 4):

| | Blank | Standard | Sample |
|-------------------------------------|-------|----------|--------|
| WR (mL) | 1,0 | 1,0 | 1,0 |
| Standard ^(Note 1,3) (μL) | -- | 10 | -- |
| Sample (μL) | -- | -- | 10 |

- Mix and incubate for 5 min. at 37°C or 10 min. at room temperature.

- Read the absorbance (A) of the samples and Standard, against the Blank. The colour is stable for at least 60 minutes.

CALCULATIONS

$$\frac{(A) \text{ Sample} - (A) \text{ Blank}}{(A) \text{ Standard} - (A) \text{ Blank}} \times 200 (\text{Standard conc.}) = \text{mg/dL cholesterol in the sample}$$

Conversion factor: mg/dL x 0,0258= mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Risk evaluation^{5,6}:

| | |
|---------------------|------------|
| Less than 200 mg/dL | Normal |
| 200-239 mg/dL | Borderline |
| ≥ 240 mg/dL | High |

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0 mg/dL to *linearity limit* of 900 mg/DL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

| | Intra-assay (n=20) | | Inter-assay (n=20) | |
|--------------|--------------------|------|--------------------|------|
| | Mean (mg/dL) | SD | CV (%) | |
| Mean (mg/dL) | 90,4 | 1,87 | 92,8 | 193 |
| SD | 1,15 | 1,01 | 1,98 | 2,39 |
| CV (%) | 1,27 | 0,54 | 2,14 | 1,24 |

Sensitivity: 1 mg/dL = 0,00152 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99541.

Regression equation: y= 0,95293x -3,020.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Hemoglobin up to 5 g/L and bilirubin up to 10 mg/dL, do not interfere^{1,2}.

A list of drugs and other interfering substances with cholesterol determination has been reported^{3,4}.

NOTES

- CHOLESTEROL CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- LCF (Lipid Clearing Factor) is integrated in the reagent.
- Calibration with the aqueous Standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- Naito H.K. Cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.
- Meiattini F. et al. The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

1. Naito H.K. Cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.
2. Meiattini F. et al. The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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| | |
|--------------|---|
| Ref: 1001090 | R1: 10 x 50 mL, R2: 10 → 50 mL, CAL: 1 x 5 mL |
| Ref: 1001091 | R1: 10 x 20 mL, R2: 10 → 20 mL, CAL: 1 x 5 mL |
| Ref: 1001092 | R1: 4 x 125 mL, R2: 4 → 125 mL, CAL: 1 x 5 mL |
| Ref: 1001093 | R1: 4 x 250 mL, R2: 4 → 250 mL, CAL: 1 x 5 mL |