**Quantitative determination of total lipids
IVD**

Store at 2-8°C

PRINCIPLE OF THE METHOD

Unsaturated lipids react with sulphuric acid to form carbonium ions. In a second step the carbonium ions react with phosphovainilline to give a pink colour.

The intensity of the color formed is proportional to the total lipids concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

The lipids are organic compounds whose more important function is the one to act like fuel.

They have an extraordinary yield, favored by the possibility of storing itself in remarkable amounts like fatty weave. Other functions: they are constituent of biological membranes, form protective fatty structures of the internal organs; provide important compounds in the formation with diverse hormones.

Great part of the interest in the study of the increase of these compounds must to the connection between hyperlipemia and arteriosclerosis, diabetes and cardiac disease^{5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Phosphovainilline Phosphoric acid	235 mmol/L 2 mmol/L
TOTAL LIPIDS CAL	Total Lipids aqueous primary standard 750 mg/dL	
Needed reagent but not provided: Sulphuric acid p.a. (H ₂ SO ₄)		

PRECAUTIONS

R: H314-Causes severe skin burns and eye damage.
Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Reagent and calibrator are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 520 nm ≥ 0,32.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 520.
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma^{1,2}:

Stability of the sample: Total lipids are stable 24 h at room temperature (15-25°C) or 3 days at 2-8°C.

PROCEDURE

- Assay conditions:
Wavelength: 520 (490-550) nm
Cuvette: 1 cm light path
Temperature: 37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a glass tubes:

	Standard	Sample
H ₂ SO ₄ (mL)	2,5	2,5
Standard ^(Note 1,2) (µL)	100	--
Sample (µL)	--	100

- Shake thoroughly using a mechanical stirrer.
- Incubate for 10 minutes in a boiling water bath (100°C).

- Cool in iced water and transfer into a cuvette:

	Blank	Standard	Sample
R (mL)	1,0	1,0	1,0
Sample Acid digest (µL)	--	--	50
Standard Acid digest (µL)	--	50	--

- Shake thoroughly using a mechanical stirrer.
- Incubate for 15 minutes at 37°C.
- Read the absorbance (A) of the samples and standard, against the Blank. The colour is stable for at least 1 hour.

CALCULATIONS

$$\frac{(A)_{\text{Sample}} - (A)_{\text{Blank}}}{(A)_{\text{Standard}} - (A)_{\text{Blank}}} \times 750 \text{ (Standard conc.)} = \text{mg/dL total lipids in the sample}$$

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum or plasma:

450 – 800 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 7,7 mg/dL to *linearity limit* of 1500 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	CV (%)	CV (%)
Mean (mg/dL)	555	919	553	919
SD	15,9	6,47	7,62	5,87
CV (%)	2,87	0,70	1,78	0,63

Sensitivity: 1 mg/dL = 0,00066 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,984

Regression equation: y=0,967x + 24,08

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

A list of drugs and other interfering substances with total lipids determination has been reported by Young et. al^{3,4}.

NOTES

- TOTAL LIPIDS CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

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