Calibrador General de Proteínas/ General Protein Calibrator

Immunoturbidimetry / Inmunoturbidimetria

Quantitative determination of Plasma Proteins

IVD
Store at 2-8°C.

PRODUCT CHARACTERISTICS

The PROT CAL is a liquid serum from human origin containing protein concentrations suitable to calibrate determinations of the constituents listed below by immunoturbidimetry method.

REAGENTS

PROT-CAL	Calibrador. Human serum. Sodium azide 0.95 g/L. Concentrations of the constituents are stated on the list below.
----------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Calibrator are traceable to Standard Reference Material ERM-DA470k/IFCC (all constituents except antithrombin III).

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Opened vials, stable for 30 days at 2-8°C.

PROCEDURE

To be used in immunoturbidimetric assays.

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LOT

This Control data sheet for lot is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number. This product was manufactured, packaged and tested in accordance to SPINREACT, S.A./S.A.U. Quality System and meets all product specifications. Our Quality System has been approved by an accredited registering body to ISO 9001.

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Immunoturbidimetry / Immunoturbidimétrie

Quantitative determination of Plasma Proteins

IVD
Store at 2-8°C.

PRODUCT CHARACTERISTICS

The PROT CAL is a liquid serum from human origin containing protein concentrations suitable to calibrate determinations of the constituents listed below by immunoturbidimetry method.

REAGENTS

PROT-CAL	Calibrator. Human serum. Sodium azide 0.95 g/L. Concentrations of the constituents are stated on the list below.
----------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Calibrator are traceable to Standard Reference Material ERM-DA470k/IFCC (all constituents except antithrombin III).

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Opened vials, stable for 30 days at 2-8°C.

PROCEDURE

To be used in immunoturbidimetric assays.

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