COVID-19 Antigen Rapid Test Kit

REF 601460

25TESTS

INTENDED USE

COVID-19 Antigen Rapid Test Kit is a solid phase immunochromatographic assay for the in vitro qualitative detection of antigen to 2019 Novel Coronavirus in human nasopharyngeal secretion or oropharyngeal secretion. This test kit provides only a preliminary test result for COVID-19 infection as a clinically-assisted diagnosis. Testing is only limited to professional laboratories .

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days. The main manifestation include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. Coronaviruses are enveloped RNA viruses that are distributed broadly among humans, other mammals, and birds and that cause respiratory, enteric, hepatic, and neurologic diseases. Seven coronavirus species are known to cause human disease. Four viruses - 229E, OC43, NL63, and HKU1 - are prevalent and typically cause common cold symptoms in immunocompetent individuals. The three other strains - severe acute respiratory syndrome coronavirus (SARS-CoV), Middle East respiratory syndrome coronavirus (MERS-CoV) and 2019 Novel Coronavirus (COVID-19) - are zoonotic in origin and have been linked to sometimes fatal illness. The COVID-19 Antigen Rapid Test Kit can detect pathogen antigens directly from nasopharyngeal swab, nasal or throat swabs specimens.

PRINCIPLE

The anti-COVID-19 Antigen Rapid Test Kit is a lateral flow immunochromatographic assay. The test uses COVID-19 antibody (test line T) and goat anti-mouse IgG (control line C) immobilized on a nitrocellulose strip. The burgundy colored conjugate pad contains colloidal gold conjugated to anti-COVID-19 antibody conjugated with colloid gold (COVID-19 conjugates) and mouse IgG-gold conjugates. When a specimen followed by assay diluent is added to the sample well, COVID-19 antigen if present, will bind to COVID-19 conjugates making antigen antibodies complex. This complex migrates through nitrocellulose membrane by capillary action. When the complex meets the line of the corresponding immobilized antibody, the complex will be combined forming a burgundy colored band which confirms a reactive test result. Absence of a colored band in the test region indicates a non-reactive test result.

The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex goat anti mouse IgG/mouse IgG-gold conjugate regardless of the color development on any of the test bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

MATERIALS SUPPLIED

Individually pouched test cassettes (25): Monoclonal anti-SARS antibodies
Single use Extraction Tubes with integrated dispensing tip (25)
Extraction Buffer (25): 380 μ L Extraction solution with detergents and reducing agents
Nylon flocked nasopharyngeal swabs (25)
Instruction for Use Sheet (1)

MATERIAL NOT PROVIDED IN KIT

Timer

Nasal swab

Viral Transport Media (Produced by Yocon or Kewei)

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use.

DO NOT FREEZE. Do not use beyond the expiration date.

WARNINGS AND PRECAUTIONS


1. For professional In Vitro diagnostic use only. Do not use after expiration date.
2. This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
3. Do not use it if the pouch is damaged or broken.
4. Test is for single use only. Do not re-use under any circumstances.
5. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
6. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
7. Humidity and temperature can adversely affect results.
8. Do not perform the test in a room with strong air flow, ie. electric fan or strong air-conditioning.
9. Test is for nasal, throat or nasopharyngeal swab specimens. If test other specimens such as serum\plasma\urine, the results are only for research.
10. Sample collection and handling procedures require specific training and guidance.
11. To obtain accurate results, use the Viral Transport Media (VTM) recommended in this Package Insert.
12. When collecting a nasal or throat swab sample, use the Nasal Swab supplied in the kit.
13. When collecting a nasopharyngeal swab sample, use a nylon flocked nasopharyngeal swab.
14. Use the appropriate Fixed Volume Pipette in accordance with test procedures.
15. Do not pour sample from the Reagent Tube into the Test Cassette sample well. Use the provided dispensing tip of Extraction Tubes or 200µL Calibrated Micropipette when adding the sample to the Test Cassette.
16. To obtain accurate results, do not use visually bloody or overly viscous samples.

SPECIMEN COLLECTION AND HANDLING

1. Nasopharyngeal Swab Sample

Use a nylon flocked nasopharyngeal swab, supplied in the kit.

To collect a nasopharyngeal swab sample, carefully insert the swab into the nostril that presents the most secretion under visual inspection. Keep the swab near the septum floor of the nose while gently pushing the swab into the posterior nasopharynx. Rotate the swab several times then remove it from the nasopharynx.

2. Nasal Swab Sample

Use the Nasal Swab, not supplied.

To collect a nasal swab sample, carefully insert the Nasal Swab (provided in the kit) into the nostril that presents the most secretion under visual inspection. Using gentle rotation, push the swab until resistance is met at the level of the turbinates (less than one inch into the nostril).

2. Oropharyngeal (throat) Swab Sample

Use the Nasal Swab, not supplied.

Insert the Nasal Swab into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, and gums. Caution has to be paid to avoid the swab to be contaminated with saliva.

SAMPLE TRANSPORT AND STORAGE

1. Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature more than 1 hour, at 2-8°C more than 8 hours. For long term storage, specimens should be kept below -20°C.
2. If transport of samples with viral transport medium (VTM) is required, minimal dilution of the sample is recommended, as dilution may result in decreased test sensitivity. Whenever possible, 1 milliliter or less is best to avoid excessive dilution of the patient sample. While holding the swab, remove the cap from the tube. Insert the swab into the tube until the breakpoint is level with the tube opening. Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. You may need to gently rotate the swab shaft to complete the breakage. Nasal, Oropharyngeal or nasopharyngeal swabs in VTM are stable for up to 72-hours at 2° to 8°C.
3. **VTM without inactivator produced by Yocon or Kewei is recommended.** Please do not use VTM with inactivator which can reduce the reactivity between antigen and antibody, other kind of VTM maybe give a false result.

TEST PROCEDURE

Allow test cassette, specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.


Note: When using viral transport medium (VTM), it is important to ensure that the VTM containing the sample is warmed to room temperature. Cold samples will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold sample to room temperature.

1. Sample preparation:

For swab sample: Prepare Extraction Tube and add all of one single buffer (380µL) into single use extraction tube. Put the swab into the Extraction Tube and spin the swab 10 times while pressing the head against the bottom and side of the Extraction Tube. Roll the swab head against the inside of the Extraction Tube as you remove it. Dispose of the used swab in biohazard waste.

For sample in viral transport medium (VTM): Prepare Extraction Tube and add 3 drops of Extraction Buffer, add 100µL sample from VTM with Micropipette.

2. Cover the extraction tube with its dispensing tip.
3. Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
4. Place the test device on a clean and level surface.
5. Add 60µL by micropipette or 2 drops using provided dispensing tip of Extraction Tubes extracted sample from the Extraction Tube to the specimen well(S) of the test device.
6. Read the results at 15 minutes. Do not interpret any result after 20 minutes.

INTERPRETATION OF RESULTS

NEGATIVE:

If only the C band is present, the absence of any burgundy color in the T band indicates

that no COVID-19 antigens are detected in the specimen. The result is negative.

POSITIVE:

In addition to the presence of C band, if T band also emerged, then the test indicates the presence of COVID-19 antigens in the specimen. The result is positive.

INVALID:

Control line fails to appear. Incorrect procedural techniques and invalid test kit are the most likely reasons for control line failure. Read the instructions carefully again and repeat the test with a new test kit. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

LIMITATIONS

1. Use fresh samples whenever possible. Frozen and thawed samples (especially repeatedly) contain particles that can block the membrane. This slows the flow of reagents and can lead to high background color, making the interpretation of results difficult.
2. Optimal assay performance requires strict adherence to the assay procedure described in this insert sheet. Deviations may lead to aberrant results.
3. A negative result for an individual subject indicates absence of detectable COVID-19 antigen. However, a negative test result does not preclude the possibility of exposure to or infection with COVID-19.
4. A negative result can occur if the quantity of the COVID-19 antigen present in the specimen is below the detection limits of the assay, or the antigen that are detected are not present during the stage of disease in which a sample is collected.
5. Positive test results do not rule out co-infections with other pathogens.
6. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis, confirmatory testing is required.

PERFORMANCE CHARACTERCS

1. Clinical Performance of Antigen Test for frozen NP swabs

The clinical performance of the device was established with a study using 182 frozen NP swabs which were originally collected in 3mL viral transport media.

		Sansure Covid-19 RT-PCR Assay		
		Positive	Negative	Total
Kewei COVID-19 Antigen Rapid Assay	Positive	71	0	71
	Negative	11	100	111
	Total	82	100	182

The assay demonstrated acceptable clinical sensitivity for frozen NP samples is 86.59% (95% CI: 77.55%-92.34%) when compared to an molecular device produced by Sansure Inc. The assay demonstrated excellent clinical specificity 100% (95% CI: 96.30%-100%). Overall Agreement was 93.96% (89.50%-96.59%).

2. Clinical Performance of Antigen Test for fresh swabs

Using 127 fresh NP swabs which were collected freshly in Extraction Buffer in the kit to evaluate the Clinical performance of the COVID-19 Antigen Rapid Test Kit. 67 Patients who presented within 7 days of symptom onset and 60 normal persons were included in the initial primary evaluation. The Result is below:

		Sansure Covid-19 RT-PCR Assay		
		Positive	Negative	Total
Kewei COVID-19 Antigen Rapid Assay	Positive	64	0	64
	Negative	3	60	63
	Total	67	60	127

The assay demonstrated acceptable clinical sensitivity for fresh swab samples is 95.52% (95% CI: 87.64%-98.47%) when compared to a molecular device produced by Sansure Inc. The assay demonstrated excellent clinical specificity 100% (95% CI: 93.98%-100%).

Overall Agreement was 97.64% (93.28%-99.19%).

3. Limit of Detection (LoD)

Limit of detection (LoD) was determined by evaluating different concentrations of heat inactivated SARS-CoV-2 virus. LoD of the SARS-CoV-2 antigen rapid kit was confirmed as 38.5 TCID₅₀/ml.

4. High-dose Hook effect

No high dose hook effect was observed when test with up to a concentration of 6.3×10⁵ TCID₅₀/ml of heat inactive SARS-CoV-2 virus from Academy of Military Sciences PLA China.

5. Cross Reactivity

Many types of pneumonia are accompanied by fever, cough and other respiratory symptoms. In order to eliminate similar clinical symptoms of other types of pneumonia effects, pneumonia mycoplasma, influenza A, pneumococcus etc. were detected for the specific assessment. Each of the organism and viruses were tested in the absence or presence of heat inactivated SARS-CoV-2 virus (308 TCID₅₀/ml), No cross-reactivity or interference was seen.

Potential Cross-Reactant	Test Concentration	Result	
		Absence of SARS-CoV-2	presence of SARS-CoV-2
Staphylococcus aureus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Streptococcus pneumoniae	1.0×10 ⁶ Cells/ml	—	+
Measles virus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Mumps virus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Adenovirus type 3	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Mycoplasma pneumoniae MP	1.0×10 ⁶ Cells/ml	—	+
Parainfluenza type 2	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Metapneumovirus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Coronavirus OC43	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Coronavirus 229E	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Bordetella parapertussis	1.0×10 ⁶ Cells/ml	—	+
Influenza B virus (Victoria series)	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Influenza B virus (Y series)	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Influenza A H1N1 (2009) virus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Influenza A H3N2 virus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Avian influenza virus H7N9	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Avian influenza virus H5N1	1.0×10 ⁵ TCID ₅₀ /ml	—	+
EB virus	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Enterovirus CA16	1.0×10 ⁵ TCID ₅₀ /ml	—	+
Rhinovirus	1.0×10 ⁵ TCID ₅₀ /ml	—	+

6. Interference

Potential Endogenous Interference of the COVID-19 Antigen Rapid Test was evaluated using natural clinical samples. No Antigen false negative or false positive results were observed with the following potential interference substances at the stated concentrations:


Human Blood	1% V/V
-------------	--------

Mucosal protein	1mg/ml
Phenylephrine	1% V/V
Oxymetazoline	1% V/V
Triamcinolone	50mg/L
Ribavirin	50mg/L

REFERENCE

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81: 85-164.
2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. Fields virology. 6th ed. Lippincott Williams & Wilkins, 2013: 825-58.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016; 24: 490-502.
4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17: 181-192.

[Symbols]

Symbol	Used for	Symbol	Used for
	Use-by date		Consult instructions for use
	Batch code		In vitro diagnostic medical device
	Temperature limit		Manufacturer
CE	CE mark		Authorized representative in the European Community
	Catalog number		Contains sufficient for<n>tests
	Date of manufacture		

BEIJING KEWEI CLINICAL DIAGNOSTIC REAGENT INC.
No.7 Yan Qi He, Xi Yi Rd., Huai Rou District, Beijing, P.R.China

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Lotus NL B.V.