**Quantitative determination of creatinine**
**IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The assay is based on the reaction of creatinine with sodium picrate as described by Jaffé.

Creatinine reacts with alkaline picrate forming a red complex. The time interval chosen for measurements avoids interferences from other serum constituents.

The intensity of the color formed is proportional to the creatinine concentration in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

Creatinine is the result of the degradation of the creatine, component of muscles, it can be transformed into ATP, that is a source of high energy for the cells. The creatinine production depends on the modification of the muscular mass, and it varies little and the levels usually are very stable.

Is excreted by the kidneys. With progressive renal insufficiency there is retention in blood of urea, creatinine and uric acid.

Elevate creatinine level may be indicative of renal insufficiency<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b> Picric Reagent	Picric acid	17,5 mmol/L
<b>R 2</b> Alkaline Reagent	Sodium hydroxide	0,29 mol/L
<b>Optional</b>	SPINTROL H CAL	

**PRECAUTIONS**

R1/ R2: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

DUAL MODE: Ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 492 nm  $\geq$  1,80.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or colorimeter measuring at 492 nm (490-510).
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

**SAMPLES**

- Serum or heparinized plasma<sup>1</sup>.
- Creatinine stability: 24 hours at 2-8°C.
- Urine (24 h)<sup>1</sup>: Dilute sample 1/50 with distilled water. Mix. Multiply results by 50 (dilution factor);
- Creatinine stability: 7 days at 2-8°C.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:

Male 0,7 - 1,4 mg/dL  $\cong$  61,8 – 123,7  $\mu$ mol/L  
 Female 0,6 - 1,1 mg/dL  $\cong$  53,0 – 97,2  $\mu$ mol/L

Urine: 15-25 mg/Kg/24 h

Male 10 - 20 mg/Kg/24 h  
 Female 8 - 18 mg/Kg/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

**APPLICATION SPINLAB 180**

Name	CREATININE	Ref. male low	0.70
Abbr. Name	CRE	Ref. male high	1.40
Mode	Twopoint	Ref. female low	0.60
Wavelength	505 nm	Ref. female high	1.10
Units	mg/dL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.20	Control 1	*
High Conc.	15.00	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	150 $\mu$ L		
rerun volume	150 $\mu$ L		
Sample			
normal volume	30.0 $\mu$ L		
rerun volume	15.0 $\mu$ L		
R2 bottle (mL)	25 mL		
normal volume	150.0 $\mu$ L		
rerun volume	150.0 $\mu$ L		
Predilución	No		
Slope blank	No		
Point one, two	24, 103 sec.		
<b>Factor</b>			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
R. Abs. Deviation	3.000 Abs		

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From *detection limit* of 0,000 mg/dL to *linearity limit* of 35 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	CV (%)	
Mean (mg/dL)	0,92	3,43	0,96	3,50
SD	0,03	0,07	0,04	0,09
CV (%)	2,76	1,90	3,97	2,51

**Sensitivity:** 1 mg/dL = 0, 0407  $\Delta$ Abs/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0, 99584

Regression equation: y= 0,953x + 0,075

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

**BIBLIOGRAPHY**

1. Murray R.L. Creatinine. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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**APPLICATION AU SPINLAB 180**

Name	CREATININE	Ref. male low	0.70
Abbr. Name	CRE	Ref. male high	1.40
Mode	Twopoint	Ref. female low	0.60
Wavelength	505 nm	Ref. female high	1.10
Units	mg/dL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.20	Control 1	*
High Conc.	15.00	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	150 µL		
rerun volume	150 µL		
Sample			
normal volume	30.0µL		
rerun volume	15.0µL		
R2 bottle (mL)	25 mL		
normal volume	150.0 µL		
rerun volume	150.0µL		
Predilución	No		
Slope blank	No		
Point one, two	24, 103 sec.		
<b>Factor</b>			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
R. Abs. Deviation	3.000 Abs		

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