### Quantitative determination of direct bilirubin IVD

Store at 2-8°C

#### PRINCIPLE OF THE METHOD

Direct bilirubin (conjugated) couples with the diazo reagent in the presence of sulfamic acid to form azobilirubin. The intensity of color formed is proportional to the bilirubin concentration in the sample tested. The increase of absorbance at 546 nm is directly proportional to the direct bilirubin concentration.

#### CLINICAL SIGNIFICANCE

Bilirubin is caused by the degradation of hemoglobin and exists in two forms. Unconjugated bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted. Hyperbilirubinemia is the result of an increase of bilirubin in plasma. Possible causes: **Total bilirubin:** Increase hemolysis, genetic, neonatal jaundice, ineffective erythropoiesis and presence of drugs. **Direct bilirubin:** Hepatic cholestasis, genetic, hepatocellular damage. Clinical diagnosis should not be made based on a single test result; it should integrate clinical and other laboratory data.

#### REAGENTS

<b>R 1</b>	Sulfamic acid	100 mM
<b>R 2</b>	2,4-DPD	0,5 mM
	Hydrochloric acid (HCl)	0,3 M

#### PRECAUTIONS

R1: H314 - Irritation or skin corrosion. / R2: H290- Corrosive to metals. H335 - May cause respiratory irritation. H314 - Irritation or skin corrosion.

R2: contains HCl and 2,4-DPD.

Follow the safety advice given in MSDS and product label.

#### PREPARATION

The reagents are provided in a ready to use format.

#### STORAGE AND STABILITY

The reagents are stable until the expiry date stated on the label when stored at

2-8°C, protected from light and contaminations are prevented during their use. Do not use reagents over the expiration date.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.

#### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 546 nm.
- General laboratory equipment.

#### SAMPLES

Serum or plasma, free of hemolysis. Protect samples from light.

Stability of the sample: 4 days at 2-8°C or 2 months at -20°C.

#### REFERENCE VALUES

Direct bilirubin 0 – 0,2 mg/dL (0 – 3,42 µmol/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

#### INTERFERENCES

Interferences from hemolysis, lipemia and ascorbic acid were evaluated for this direct bilirubin method on a Spintech 240 analyzer. Two concentrations of direct bilirubin were evaluated. No interferences were observed for lipemia (Intralipid) up to 350 mg/dL and ascorbic acid up to 40 mg/L. Hemolysis causes decreased direct bilirubin values.

A list of drugs and other interfering substances with bilirubin has been reported by Young et. al <sup>4,5</sup>.

#### APPLICATION SPINLAB 180

Name	BILIRUBIN D	Ref. male low	*
Abbr.name	BILD	Ref. male high	*
Mode	Twopoint	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/dL	Ref. Ped. low	*
Decimals	2	Ref. Ped. high	*
Low conc.	0,03 mg/dL	Valor pánico bajo	*
High conc.	9,00 mg/dL	Valor pánico alto	*
Calibrator	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
Normal volume	240 µL	Normal volume	** µL
Rerun volume	240 µL	Rerun volume	** µL
Sample		Sample	
Normal volume	15,0 µL	Normal volume	** µL
Rerun volume	15,0 µL	Rerun volume	** µL
R2 bottle (mL)	5 mL		
Normal volume	60 µL		
Rerun volume	60 µL		
Predilution	No		
Incubation time	-3.236 seg.	Incubation time	**
Factor	**	Factor	**
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs Abs	Low Absorbance	-0.100 Abs Abs
High Absorbance	3.000 Abs Abs	High Absorbance	3.000 Abs Abs
R.Abs.L.Limit	-0.100 Abs Abs	R.Abs.L.Limit	-0.100 Abs Abs
R.Abs.H.Limit	3.000 Abs Abs	R.Abs.H.Limit	3.000 Abs Abs
R. Abs. Deviation	3.000 Abs Abs	R. Abs. Deviation	3.000 Abs Abs

#### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 0,03 mg/dL to *linearity limit* of 9 mg/dL. If the results obtained are greater than the linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

Mean (mg/dL)	Inter assay (n= 40)		Intra assay (n= 80)	
	0,7458	2,444	0,7458	2,444
SD	0,05868	0,0550	0,0276	0,024
CV (%)	7,9	2,2	3,7	1,0

**Sensitivity:** 1 mg/dL = 0,040 Abs. units

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x) on a Spintech 240 analyzer. The results obtained using 53 samples ranging from 0,06 a 9 mg/dL (1,02 to 153,9 µmol/L) were:

Correlation coefficient (r): 0,9986

Regression equation:  $y = 1,0056x - 0,1046$

The results of the performance characteristics depend on the analyzer used.

#### NOTES

1. SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

#### BIBLIOGRAPHY

1. David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307–328.
2. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2: 481-491.
3. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307–328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2: 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP1001047	Cont.	R 1: 10 x 20 mL
		R 2: 10 x 5 mL

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**APPLICATION AU SPINLAB 180**

Name	BILIRUBIN D	Ref. male low	*
Abbr.name	BILD	Ref. male high	*
Mode	Twopoint	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/dL	Ref. Ped. low	*
Decimals	2	Ref. Ped. high	*
Low conc.	0.03 mg/dL	Valor pánico bajo	*
High conc.	9.00 mg/dL	Valor pánico alto	*
Calibrator	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
Normal volume	240 µL	Normal volume	** µL
Rerun volume	240 µL	Rerun volume	** µL
Sample		Sample	
Normal volume	15.0 µL	Normal volume	** µL
Rerun volume	15.0 µL	Rerun volume	** µL
R2 bottle (mL)	5 mL		
Normal volume	60 µL		
Rerun volume	60 µL		
Predilution	No		
Incubation time	-3.236 seg.	Incubation time	**
Factor	**	Factor	**
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs Abs	Low Absorbance	-0.100 Abs Abs
High Absorbance	3.000 Abs Abs	High Absorbance	3.000 Abs Abs
R.Abs.L.Limit	-0.100 Abs Abs	R.Abs.L.Limit	-0.100 Abs Abs
R.Abs.H.Limit	3.000 Abs Abs	R.Abs.H.Limit	3.000 Abs Abs
R. Abs. Deviation	3.000 Abs Abs	R. Abs. Deviation	3.000 Abs Abs

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**BIBLIOGRAPHIE**

1. David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
2. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
3. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Text book of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP1001047	Cont.	R 1: 10 x 20 mL
		R 2: 10 x 5 mL

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### BIBLIOGRAFIA

- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP1001047	Cont.	R 1: 10 x 20 mL
		R 2: 10 x 5 mL