Lp(a)-turbilatex

Latex turbidimetry

Quantitative determination of Lipoprotein (a) (Lp(a)) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The Lp(a)-turbilatex is a quantitative turbidimetric test for the measurement of Lp(a) in human serum or plasma. Latex particles coated with antibodies anti-Lp(a) are agglutinated when mixed with samples containing Lp(a). The agglutination causes an absorbance change, dependent upon the Lp(a) contents of sample that can be quantified by comparison from a calibrator of known Lp(a) concentration.

CLINICAL SIGNIFICANCE

Lp(a) is a low density lipoprotein-like particle containing apolipoprotein B-100 disulphide-linked to one large glycoprotein called apolipoprotein (a). Many investigators have confirmed that a high Lp(a) concentration represents an indicator of risk for cardiovascular disease, especially when serum LDL-cholesterol or Apo B are elevated. The quantification of Lp(a) in serum or plasma is important for identification of individuals at risk for developing atherosclerosis.

REAGENTS

Diluent (R1)	Glycine buffer 50 mmol/L. Preservative.
Latex (R2)	Latex particles coated with rabbit polyclonal anti-human Lp(a). Preservative.
Optional	Ref.: 1107022 Lp(a) Calibrator. Ref.: 1107024 Lp(a) Control.

CALIBRATION

The sensitivity of the assay and the target value of the calibrator have been standardized against an Internal Reference Material. It is not recommended the use of other commercially available Lp(a) calibrators.

PREPARATION

Calibration Curve Prepare the following Lp(a) calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the Lp(a) calibrator by the corresponding factor stated in table below to obtain the Lp(a) concentration of each dilution.

Calibrator dilution	1	2	3	4	5
Lp(a) Calibrator (µL)	--	25	50	75	100
NaCl 9 g/L (µL)	100	75	50	25	-
Factor	0	0,25	0,5	0,75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date. Do not freeze; frozen Latex or Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Normal values up to 30 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Lp(a) Control Ref.: 1107024.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINLAB 180 APPLICATION

Name	Lp(a)	Ref. male low	0.0 mg/dL
Abbr. Name	Lp(a)	Ref. male high	30.0 mg/dL
Mode	Twopoints	Ref. female low	0.0 mg/dL
Wavelength	578 nm	Ref. female high	30.0 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	12.0 mg/dL	Control 1	*
High Conc.	80.0 mg/dL	Control 2	*
Calibrator name	CAL Lp(a)	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µ L		
rerun volume	240 µ L		
Sample			
normal volume	5.0 µL		
rerun volume	3.0 µ L		
R2 bottle (mL)	5 mL		
normal volume	60 µ L		
rerun volume	60 µ L		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Reagent blank	Yes		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

INTERFERENCES

Hemoglobin (1000 mg/dL), bilirubin (60 mg/dL) and lipemia (10 g/L) do not interfere. Other substances may interfere⁵.

PERFORMANCE CHARACTERISTICS

- Linearity:** Up to 100 mg/dL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depends on the sample to reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Quantification limit:** Values less than 3,17 mg/dL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 650 mg/dL.
- Precision:**

	Level 1	Level 2	Level 3
Mean (mg/dL)	18,81	34,12	47,60
SD	0,59	0,48	0,33
CV	3,15	1,42	0,69

- Accuracy:** The Spinreact method was compared with another manufacturer. The study was performed with 35 serum samples. Both tests were performed on a Spintech240. Both tests were calibrated with their respective calibrators. The correlation coefficient (r) was 0,9815 and the regression equation $y = 1,5734x + 1,5873$. The results of the performance characteristics depend on the analyzer used.

NOTES

- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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APPLICATION AU SPINLAB 180

Name	Lp(a)	Ref. male low	0.0 mg/dL
Abbr. Name	Lp(a)	Ref. male high	30.0 mg/dL
Mode	Twopoints	Ref. female low	0.0 mg/dL
Wavelength	578 nm	Ref. female high	30.0 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	12.0 mg/dL	Control 1	*
High Conc.	80.0 mg/dL	Control 2	*
Calibrator name	CAL Lp(a)	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	5.0 µL		
rerun volume	3.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Reagent blank	Yes		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

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