Quantitative determination of immunoglobulin E (IgE)
IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The IgE-Turbilatex is a quantitative turbidimetric test for the measurement of IgE in human serum or plasma.

Latex particles coated with mouse IgG anti-human IgE are agglutinated when mixed with samples containing IgE. The agglutination causes an absorbance change, dependent upon the IgE contents of the patient sample that can be quantified by comparison from a calibrator of known IgE concentration.

CLINICAL SIGNIFICANCE

IgE is an immunoglobulin with a molecular weight of approximately 190,000 normally present in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment assessment of disease progression, or postoperative prognosis for such conditions.

REAGENTS

Diluent (R1)	Glycine buffer, pH 8.3. Preservative.
Latex (R2)	Latex particles coated with mouse IgG anti-human IgE, pH 7.3. Preservative.
Optional	Ref: 1107052 IgE CAL (Calibrator) Ref: 1107053 IgE Control.

CALIBRATION

The sensitivity of the assay and the target value of the calibrator have been standardized against the International Reference Preparation of IgE 75/502 (2nd IRP 1981) from WHO. It is not recommended the use of other commercially available IgE calibrators.

PREPARATION

Reagents: Ready to use.

Calibration curve: Prepare the following IgE calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the IgE calibrator by the corresponding factor stated in table below to obtain the IgE concentration of each dilution.

Calibrator dilution	1	2	3	4	5
IgE Calibrator (µL)	--	12.5	25	50	100
NaCl 9 g/L (µL)	100	87.5	75	50	--
Factor	0	0.125	0.25	0.5	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Frozen Latex and Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Up to 350 IU/mL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact Control Serum IgE is available (Ref.: 1107053).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should be not made on findings of a single test result but should integrate both clinical and laboratory data.

SPINLAB 180 APPLICATION

Name	IgE	Ref. male low	0 UI/mL
Abbr. Name	IgE	Ref. male high	300 UI/mL
Mode	Twopoints	Ref. female low	0 UI/mL
Wavelength	578 nm	Ref. female high	300 UI/mL
Units	UI/mL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	25 UI/mL	Control 1	*
High Conc.	1000 UI/mL	Control 2	*
Calibrator name	CAL IgE	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		CALIBRATION	
Sample blank	No	Calibrator	CAL IgE
R1 bottle (mL)	25 mL	Repeat	2
normal volume	217 µL	Interval	0
rerun volume	217 µL	Max. Inaccuracy	0.01%
Sample			
normal volume	5.0 µL		
rerun volume	3.0 µL		
R2 bottle (mL)	5 mL		
normal volume	116 µL		
rerun volume	116 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Factor	*		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 1000 IU/mL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depend on the sample to reagent / ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Limit detection:** Values less than 3,04 IU/mL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 60,000 IU/mL.
- Sensitivity:** Δ 0.5 mA. IU/mL.
- Precision:**

Within run			
n= 21	Mean (IU/mL)	SD (IU/mL)	CV (%)
	64,0	1,07	1,68
	121,0	1,43	1,18
	295,4	1,20	0,41
	491,3	4,01	0,82
	739,8	2,63	0,36

Between days			
n= 21	Mean (IU/mL)	SD (IU/mL)	CV (%)
	61,0	2,9	4,76
	406	7,8	1,92

- Accuracy:** Results obtained using this reagents (y) were compared to those obtained using a nephelometric method (x). 92 samples were assayed by both methods. The correlation coefficient (r) was 0,9748 and the regression line equation $y = 1,0321x + 8,7369$.

The results of the performance characteristics depend on the analyzer used.

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APPLICATION AU SPINLAB 180

Name	IgE	Ref. male low	0 UI/mL
Abbr. Name	IgE	Ref. male high	300 UI/mL
Mode	Twopoints	Ref. female low	0 UI/mL
Wavelength	578 nm	Ref. female high	300 UI/mL
Units	UI/mL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	25 UI/mL	Control 1	*
High Conc.	1000 UI/mL	Control 2	*
Calibrator name	CAL IgE	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		CALIBRATION	
Sample blank	No	Calibrator	CAL IgE
R1 bottle (mL)	25 mL	Repeat	2
normal volume	217 µL	Interval	0
rerun volume	217 µL	Max. Inaccuracy	0.01%
Sample			
normal volume	5.0 µL		
rerun volume	3.0 µL		
R2 bottle (mL)	5 mL		
normal volume	116 µL		
rerun volume	116 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Factor	*		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		

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Name	IgE	Ref. male low	0 UI/mL
Abbr. Name	IgE	Ref. male high	300 UI/mL
Mode	Twopoints	Ref. female low	0 UI/mL
Wavelength	578 nm	Ref. female high	300 UI/mL
Units	UI/mL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	25 UI/mL	Control 1	*
High Conc.	1000 UI/mL	Control 2	*
Calibrator name	CAL IgE	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		CALIBRATION	
Sample blank	No	Calibrator	CAL IgE
R1 bottle (mL)	25 mL	Repeat	2
normal volume	217 µL	Interval	0
rerun volume	217 µL	Max. Inaccuracy	0.01%
Sample			
normal volume	5.0 µL		
5.0 µL			
rerun volume	3.0 µL		
3.0 µL			
R2 bottle (mL)	5 mL		
normal volume	116 µL		
rerun volume	116 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Factor	*		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		

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