**Quantitative determination of anti-streptolysin O (ASO)**
**IVD**

Store 2-8°C.

**PRINCIPLE OF THE METHOD**

The ASO-Turbilatex is a quantitative turbidimetric test for the measurement of ASO in human serum or plasma.

Latex particles coated with streptolysin O (SLO) are agglutinated when mixed with samples containing ASO. The agglutination causes an absorbance change, dependent upon the ASO contents of the patient sample that can be quantified by comparison from a calibrator of known ASO concentration.

**CLINICAL SIGNIFICANCE**

SLO is a toxic immunogenic exoenzyme produced by  $\beta$ -hemolytic Streptococci of groups A, C and G. Measuring the ASO antibodies are useful for the diagnostic of rheumatoid fever, acute glomerulonephritis and streptococcal infections. Rheumatic fever is an inflammatory disease affecting connective tissue from several parts of human body as skin, heart, joints etc... and acute glomerulonephritis is a renal infection that affects mainly to renal glomerulus.

**REAGENTS**

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8.2. Preservative.
<b>Latex (R2)</b>	Latex particles coated with streptolysin O, pH 10.0. Preservative.
<b>ASO-CAL</b>	Calibrator. Human serum. ASO concentration is stated on the vial label.
<b>Optional</b>	Ref.:1102114 Control ASO/CRP/RF Level L Ref.:1102115 Control ASO/CRP/RF Level H

**PRECAUTIONS**

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

**CALIBRATION**

Use ASO Calibrator Reference 1107012.

The sensitivity of the assay and the target value of the calibrator have been standardized against the ASO International Calibrator from NIBSC ASO. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

**PREPARATION**

Ready for use.

**ASO Calibrator:** Reconstitute (→) with 1.0 mL of distilled water. Mix gently and incubate at room temperature for 10 minutes before use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** Presence of particles and turbidity.

**ASO Calibrator:** Stable for 1 month at 2-8°C or 3 months at -20°C.

**Do not freeze; frozen Latex or Diluent could change the functionality of the test.**

**ADDITIONAL EQUIPMENT**

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

**SAMPLES**

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Controls ASO/CRP/RF Level L (Ref.: 1102114) and Level H (Ref.: 1102115).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES**

Normal values up to 200 IU/mL (adults) and 100 IU/mL (children < 5 years old)<sup>6</sup>.

Each laboratory should establish its own reference range.

**SPINLAB 180 APPLICATION**

Name	ASO	Ref. male low	*
Abbr. Name	ASO	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	IU/mL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	20 IU/mL	Panic value low	*
High Conc.	800 IU/mL	Panic value high	*
Calibrator name	CAL ASO	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000

<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 $\mu$ L		
rerun volume	240 $\mu$ L		
Sample			
normal volume	3.0 $\mu$ L		
rerun volume	2.0 $\mu$ L		
R2 bottle (mL)	5 mL		
normal volume	60 $\mu$ L		
rerun volume	60 $\mu$ L		
Predilución	No		
Slope blank	No		
Point one,two	6, 130 sec.		

Reagent blank	Yes
Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs
Substr.Depletion	3.000 Abs

**INTERFERENCES**

Bilirubin (20 mg/dL), hemoglobin (10 g/L), lipemia (10 g/L) and rheumatoid factors (600 IU/mL), do not interfere. Other substances may interfere<sup>7</sup>.

**PERFORMANCE CHARACTERISTICS**

- Linearity limit:** Up to 800 IU/mL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/3 in NaCl 9 g/L and retested again. The linearity limit depends on the sample-reagent ratio, as well the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 20 IU/mL give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 1000 IU/mL.
- Sensitivity:**  $\Delta$  0.73 mA. IU/mL.
- Precision:** The reagent has been tested for 20 days, using three different ASO concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 100 IU/mL	+/- 200 IU/mL	+/- 400 IU/mL
Total	6.4%	5.7%	5.1%
Within Run	2.4%	1.7%	1.4%
Between Run	3.6%	4.2%	4.9%
Between Day	4.7%	3.5%	0.7%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 60 samples of different concentrations of ASO were assayed. The correlation coefficient (r) was 0.99 and the regression equation  $y = 0.915x - 4.844$ . The results of the performance characteristics depend on the analyzer used.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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**APPLICATION AU SPINLAB 180**

Name	ASO	Ref. male low	*
Abbr. Name	ASO	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	IU/mL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	20 IU/mL	Panic value low	*
High Conc.	800 IU/mL	Panic value high	*
Calibrator name	CAL ASO	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	3.0 µL		
5.0 µL			
rerun volume	2.0 µL		
3.0 µL			
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 130 sec.		
Reagent blank	Yes		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		

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