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- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarrera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
- Wibell L et al. Nephron 1973; 10: 320-331.
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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Ref.: SP1107030

Cont.

R1. Diluyente: 2 x 20 mL  
 R2. Látex: 2 x 5 mL  
 β<sub>2</sub>-m CAL: 1 x 1 mL

β<sub>2</sub>-m TURBI

# β<sub>2</sub>-m turbilatex

Latex turbidimetry

## Quantitative determination of β<sub>2</sub>-microglobulin (β<sub>2</sub>-m)

### IVD

Store 2 - 8°C.

### PRINCIPLE OF THE METHOD

The β<sub>2</sub>-m Turbilatex is a quantitative turbidimetric test for the measurement of β<sub>2</sub>-microglobulin (β<sub>2</sub>-m) in human serum, plasma or urine. Latex particles coated with anti-human β<sub>2</sub>-m are agglutinated when mixed with samples containing β<sub>2</sub>-m. The agglutination causes an absorbance change, dependent upon the β<sub>2</sub>-m contents of the patient sample that can be quantified by comparison from a calibrator of known concentration.

### CLINICAL SIGNIFICANCE

β<sub>2</sub>-m is a protein located on the surface of human lymphocytes and other nucleated cells. Free β<sub>2</sub>-m is filtered by the glomerulus and subsequently reabsorbed in the proximal tubular cells. Increased urinary excretion of β<sub>2</sub>-m is a sensitive indicator of renal insufficiency. Also, the β<sub>2</sub>-m level in serum is a useful marker of other diseases including carcinomas, lymphoid tumors, rheumatoid arthritis and AIDS.

### REAGENTS

|                                     |  |
|-------------------------------------|--|
| <b>β<sub>2</sub>-m Diluent (R1)</b> | Tris buffer 20 mmol/L, pH 8.2. Preservative.                                       |
| <b>β<sub>2</sub>-m Latex (R2)</b>   | Particles coated with goat IgG anti-human β <sub>2</sub> -m, pH 8.0. Preservative. |
| <b>β<sub>2</sub>-m CAL</b>          | Calibrator. β <sub>2</sub> -m concentration is stated on the vial label.           |
| <b>Optional:</b>                    | Ref: 1107034 β <sub>2</sub> -m Control.  |

### PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

### CALIBRATION

Use β<sub>2</sub>- Microglobulin Calibrator Reference 1107032. The sensitivity of the assay and the target value of the calibrator have been standardized against the 1<sup>st</sup> International ERM-DA470k/IFCC. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

### PREPARATION

Ready for use.

#### β<sub>2</sub>-m Calibrator:

Serum method: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.  
Urine method: Dilute reconstituted calibrator 1/6 with NaCl 9 g/L (50 μL calibrator + 250 μL NaCl 9 g/L).

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles and turbidity.

**β<sub>2</sub>-m Calibrator:** Stable for 1 month at 2-8°C or 3 months at -20°C.

### ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

### SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Fresh urine. Adjust samples to pH 7-8 by the addition of K<sub>2</sub>HPO<sub>4</sub>. Stable 2 days at 2-8°C or 2 months at -20°C.

The samples with particles or fibrin should be centrifuged before testing. Do not use hemolyzed or lipemic samples.

### QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT β<sub>2</sub>-m Control (Ref: 1107034). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

**Serum:** from 1.0 to 3.0 mg/L.

**Urine:** from 0.1 to 0.3 mg/L.

Each laboratory should establish its own reference range

## SPINLAB 180 APPLICATION

|                 |            |                  |       |
|-----------------|------------|------------------|-------|
| Name            | b2M        | Ref. male low    | *     |
| Abbr. Name      | b2M        | Ref. male high   | *     |
| Mode            | Twopoints  | Ref. female low  | *     |
| Wavelength      | 546 nm     | Ref. female high | *     |
| Units           | mg/L       | Ref. Ped. Low    | *     |
| Decimals        | 2          | Ref. Ped. High   | *     |
| Low Conc.       | 0.22 mg/L  | Panic value low  | *     |
| High Conc.      | 18.00 mg/L | Panic value high | *     |
| Calibrator name | CAL b2M    | Control 1        | *     |
| Prozone check   | No         | Control 2        | *     |
|                 |            | Control 3        | *     |
|                 |            | Correlat. factor | 1.000 |
|                 |            | Correlat. offset | 0.000 |

|                  |             |  |  |
|------------------|-------------|--|--|
| <b>DUAL MODE</b> |             |  |  |
| Sample blank     | No          |  |  |
| R1 bottle (mL)   | 25 mL       |  |  |
| normal volume    | 240 μL      |  |  |
| rerun volume     | 240 μL      |  |  |
| Sample           |             |  |  |
| normal volume    | 3.0 μL      |  |  |
| rerun volume     | 3.0 μL      |  |  |
| R2 bottle (mL)   | 5 mL        |  |  |
| normal volume    | 60 μL       |  |  |
| rerun volume     | 60 μL       |  |  |
| Predilución      | No          |  |  |
| Slope blank      | No          |  |  |
| Point one,two    | 6, 183 sec. |  |  |

|                  |            |
|------------------|------------|
| Reagent blank    | No         |
| Low Absorbance   | -0.100 Abs |
| High Absorbance  | 3.000 Abs  |
| R. Abs. L. Limit | -0.100 Abs |
| R. Abs. H. Limit | 3.000 Abs  |
| Substr.Depletion | 3.000 Abs  |

### INTERFERENCES

**Serum method:** bilirubin (20 mg/L), hemoglobin (10 g/L) and lipids (10 g/L), do not interfere. Rheumatoid factors (150 IU/mL), interfere.

**Urine method:** urea (urine) (50 g/L), uric ac. (20 g/L) and glucose (100 g/L), do not interfere.

Other substances may interfere<sup>7</sup>.

### PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 18 mg/L (serum) and 3 mg/L (urine), under the described assay conditions. Samples higher results should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity depends on the sample-reagent ratio, as well as the analyzer used. It will be higher by decreasing sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0.09 mg/L (serum) and 0.01 mg/L (urine) give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 100 mg/L (serum) and 20 mg/L (urine).
- Sensitivity:** Δ 0.048 A. mg/L (serum) and Δ 0.228 A. mg/L (urine).
- Precision:** The reagent has been tested for 20 days, using three different β<sub>2</sub>-m concentrations in a EP5-based study.

| EP5         | CV (%)     |              |              |
|-------------|------------|--------------|--------------|
|             | +/- 1 mg/L | +/- 3.2 mg/L | +/- 8.5 mg/L |
| Total       | 4.0%       | 3.4%         | 1.7%         |
| Within Run  | 2.8%       | 2.0%         | 1.2%         |
| Between Run | 1.7%       | 1.5%         | 1.2%         |
| Between Day | 2.2%       | 2.4%         | 0.0%         |

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 36 samples of different concentrations of β<sub>2</sub>-m were assayed. The correlation coefficient (r) was 0.97 and the regression equation y = 1.709x - 2.627. The results of the performance characteristics depend on the analyzer used.

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|                  |             |                  |       |
|------------------|-------------|------------------|-------|
| Name             | b2M         | Ref. male low    | *     |
| Abbr. Name       | b2M         | Ref. male high   | *     |
| Mode             | Twopoints   | Ref. female low  | *     |
| Wavelength       | 546 nm      | Ref. female high | *     |
| Units            | mg/L        | Ref. Ped. Low    | *     |
| Decimals         | 2           | Ref. Ped. High   | *     |
| Low Conc.        | 0.22 mg/L   | Panic value low  | *     |
| High Conc.       | 18.00 mg/L  | Panic value high | *     |
| Calibrator name  | CAL b2M     | Control 1        | *     |
| Prozone check    | No          | Control 2        | *     |
|                  |             | Control 3        | *     |
|                  |             | Correlat. factor | 1.000 |
|                  |             | Correlat. offset | 0.000 |
| <b>DUAL MODE</b> |             |                  |       |
| Sample blank     | No          |                  |       |
| R1 bottle (mL)   | 25 mL       |                  |       |
| normal volume    | 240 µL      |                  |       |
| rerun volume     | 240 µL      |                  |       |
| Sample           |             |                  |       |
| normal volume    | 3.0 µL      |                  |       |
| rerun volume     | 3.0 µL      |                  |       |
| R2 bottle (mL)   | 5 mL        |                  |       |
| normal volume    | 60 µL       |                  |       |
| rerun volume     | 60 µL       |                  |       |
| Predilución      | No          |                  |       |
| Slope blank      | No          |                  |       |
| Point one,two    | 6, 183 sec. |                  |       |
| Reagent blank    | No          |                  |       |
| Low Absorbance   | -0.100 Abs  |                  |       |
| High Absorbance  | 3.000 Abs   |                  |       |
| R. Abs. L. Limit | -0.100 Abs  |                  |       |
| R. Abs. H. Limit | 3.000 Abs   |                  |       |
| Substr.Depletion | 3.000 Abs   |                  |       |

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**BIBLIOGRAFIA**

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarrera M et al. Digestive Diseases and Sciences 1997; 42: 762-766
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93
- Wibell L et al. Nephron 1973; 10: 320-331
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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|                 |                                 |
|-----------------|---------------------------------|
| Ref.: SP1107030 | R1 Diluent :2 x 20 mL           |
|                 | R2 Latex :2 x 5 mL              |
|                 | β <sub>2</sub> -m CAL :1 x 1 mL |

Cont.